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









INTENDED USE

- The bracket and disposable guide provide physicians with a tool for performing needle-guided (or catheter) procedures with the use of diagnostic ultrasound transducers.

WARNING

- In the U.S.A., federal law restricts device to sale or use by or on order of a physician.
- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's manual.

I. ATTACHING BRACKET TO TRANSDUCER _____

WARNING

- Bracket is packaged non-sterile and is reusable. To avoid possible patient contamination, ensure bracket is properly cleaned, sterilized or disinfected before each use.
- If bracket does not fit properly, do not use.

1. Attach bracket to transducer by aligning locating features. Ensure bracket is firmly attached.

NOTE

II. SELECTING BRACKET ANGLE _____

- For multi-angle brackets only.
1. Pull up lock pin.
 2. Move to desired angle.
 3. Press lock pin securely into position.

WARNING

III. COVERING THE TRANSDUCER AND BRACKET _____

- Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.
- Disposable components are packaged sterile and are single-use only. Do not use if integrity of packaging is violated or if expiration date has passed.
- For illustration purposes only, transducer may be shown without a transducer cover. Always place a cover over transducer to protect patients and users from cross-contamination.

1. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
2. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
3. Secure cover with enclosed bands.
4. Inspect cover to ensure there are no holes or tears.

IV. ATTACHING NEEDLE GUIDE TO BRACKET _____

NOTE

- L006A and V013A: If using the 53°/55° or 64° angles, do not set the lock pin until after the disposable needle guide has been attached. Attach the needle guide at the 42°/46° angle, then select the desired angle, 53°/55° or 64°.

1. Using proper sterile technique, snap unlocked needle guide onto attachment area of bracket.
2. Push lock into locked position.
3. Select appropriate size needle guide insert and slide into position.

NEEDLE PATH VERIFICATION _____

4. Before using component for first time, perform needle path verification to verify system and biopsy guide relationships, as described in system manual.
5. For multi-angle brackets, verify bracket angle corresponds with system software before each use.
6. Use appropriate needle length to reach target area.

V. ACTIVATING NEEDLE GUIDE QUICK-RELEASE

1. Press tab toward bracket to activate needle guide quick-release function.

VI. DISASSEMBLING BRACKET FOR CLEANING

NOTE | • See illustrations on page 8.

VII. DISPOSAL

WARNING | • Dispose of single-use components as infectious waste. Clean and sterilize reusable components after each use.

VIII. CLEANING, DISINFECTING AND STERILIZING BRACKET

WARNING | • Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

- See your system’s manual for reprocessing transducer between use.
- The recommended protocol for reprocessing the reusable bracket has been evaluated for effectiveness and compatibility. If choosing to use a reprocessing protocol other than those listed, a validation must first be completed to ensure effectiveness and compatibility.
- Do not gas sterilize or autoclave bracket.

CLEANING THE BRACKET

1. After each use, remove bracket and needle guide from transducer, discarding single-use, disposable components. Remove visible contaminants from bracket surface using a small, soft instrument brush. Keep bracket from drying out until complete cleaning can be accomplished.
2. Soak bracket for minimum of five minutes in neutral pH, low foaming enzymatic detergent, such as ENZOL® Enzymatic Detergent (Johnson & Johnson). While immersed, use instrument brush to remove trapped contaminants from surfaces and features. If visible contaminants cannot be easily removed, repeat soaking procedure for an additional five minutes. Remove bracket from cleaning solution and remove any remaining residue with dry wipe. Follow cleaning solution manufacturer’s directions for use and recommendations for concentration.

DISINFECTING AND STERILIZING BRACKET

- High-level disinfect bracket using CIDEX OPA® ortho-Phthalaldehyde Solution (Johnson & Johnson) or equivalent .55% ortho-phthalaldehyde-based solution. Follow manufacturer’s instructions and recommendations for concentration, time of contact and post-process procedure.
- High-level disinfect or sterilize bracket using CIDEX® Activated Dialdehyde Solution (Johnson & Johnson) or equivalent 2% glutaraldehyde-based solution, or a hydrogen peroxide-based solution. Follow manufacturer’s instructions and recommendations for concentration, time of contact and post-process procedure.
- Inspect bracket for damage such as cracks or breakage. If damage is evident, discontinue use of bracket and contact your ultrasound system support representative.

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COPYRIGHT © 2006 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ULTRA-PRO II IS A TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. U.S. PATENTS NOS. 5,941,889 AND 6,361,499.
MPS GmbH