
Percutaneous Instrument Guide


For use with CIVCO Classic and
Multi-Purpose Workstation™ Steppers


Reference Guide

Percutaneous Instrument Guide


Mode D -Two Plates
D-funktion – to plader
Modus D -twee platen
D-tila - kaksi levyä
Mode D – Deux plaques
Modus D - zwei Platten
Θέση D - Δύο πλάκες
Modalità D - due piastre
モードD - 2プレート
Modo D -Duas placas
Modo D -Dos placas
Läge D -Två plattor


Mode F -One Plate
F-funktion – én plade
Modus F -één plaat
F-tila - yksi levy
Mode F – Une plaque
Modus F - eine Platte
Θέση F - Μία πλάκα
Modalità F - una piastra
モードF - 1プレート
Modo F -Uma placa
Modo F -Una placa
Läge F -En platta

Percutaneous Instrument Guide

- I. Removing plate
Fjernelse af pladen
Plaat verwijderen
Levyn irrottaminen
Retrait de la plaque
Entfernen der Platte
Αφαίρεση πλάκας
Rimozione piastra
プレートの取り外し
Remover a placa
Extracción de la placa
Borttagning av platta

- II. Mounting guide to platform
Montering af guiden på platformen
Montering af guiden på platformen
De voerder aan het platform bevestigen
Ohjaimen kiinnittäminen alustaan
Montage du guide sur la plate-forme
Befestigen der Führung an der Plattform
Τοποθέτηση του οδηγού στη βάση
Fissaggio della guida alla piattaforma
プレートへのガイドの取り付け
Montar a guia na plataforma
Montaje de la guía en la plataforma
Montering av guiden på plattformen

INTENDED USE | • *This device is intended to be used to aid in the percutaneous positioning and placement of needles or needlelike therapeutic probes.*

WARNING | • *In the USA, federal law restricts device to sale or use by or on order of a physician.*
• *Before use, you should be trained in cryotherapy. For instructions on preparation, safe operation and care of your transducer and cryotherapy system, see your system’s manual.*

NOTE | • Disposable components are packaged sterile and are single-use only. Do not use if integrity of packaging is violated or if expiration date has passed. Percutaneous Instrument Guide (guide) cannot be resterilized after use.
• Do not use guide if it appears to be damaged.
• Ensure guide is fully engaged and fits securely to the platform.
• Guide to be used with 2.4mm therapeutic probe.
• Only to be used with compatible CIVCO Steppers.

ASSEMBLING TEMPLATE _____

1. Using sterile technique, remove guide from package.
2. Configure guide for desired mode, one plate or two (*Refer to illustrations*).
3. Disengage clip to remove plate for Mode F only (*Refer to illustration I*).
4. Mount guide to platform (*Refer to illustration II*).
5. Ensure guide is secured to stepper platform (*Refer to illustration II*).

DISPOSAL _____

NOTE | • Dispose of single-use components as infectious waste.

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COPYRIGHT © 2006 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. U.S. PRINTED IN USA.