
# Accessories Kit

*Drape, Grid, NeoGuard™ Covers*

**Reference Guide**



CE  
0120

	<p>Storage conditions Opbevaring Opslagvoorwaarden Varastointioloisheet</p>	<p>Conditions de conservation Bedingungen zur Aufbewahrung Όοί έπεάο όπεάί ςο Indicazioni per la conservazione</p>	<p>保管条件 Condições de armazenamento Condiciones de almacenamiento Förvaringsförhållanden</p>
---	---	--	---

	<p>UV light sensitive UV lysfølsom Uv lichtsensitief Ultraviolettalolle herkka</p>	<p>Sensible aux rayons UV UV-lichtempfindlich αόαβοεςοί όοί οόάνεπαάο οοο Sensibile ai raggi ultravioletti</p>	<p>紫外線に感光し止す Sensível à luz UV Sensible a rayos ultravioleta Känslig mot UV-ljus</p>
---	--	--	--




- INTENDED USE** | • Provides physicians with a series of precisely placed needle guide channels used to facilitate accurate placement and insertion in relation to the ultrasound imaging probe.
- WARNING** | • In the USA, federal law restricts device to sale or use by or on order of a physician.  
 • Before use, you should be trained in brachytherapy. For instructions on preparation, safe operation and care of your transducer and brachytherapy system, see your system's manual.  
 • Do not use if integrity of packaging is violated or if expiration date has passed.  
 • Do not use grid if it appears to be cracked or broken.  
 • Grid cannot be resterilized after use.  
 • Ensure graphics on ultrasound monitor match graphics on grid.  
 • Take care when inserting needles into grid to avoid breaking needle.

**ASSEMBLING GRID**

1. Using sterile technique, remove grid from package.
2. Select appropriate graphics on grid to match ultrasound monitor.
3. Tighten locking knobs.

**DISPOSAL**

- WARNING** | • Dispose of single-use components as infectious waste.

**INTENDED USE**

- The cover allows use of the transducer in scanning and needle guided procedures for body surface, endocavity and intra-operative diagnostic ultrasound, while helping to prevent transfer of microorganisms, body fluids and particulate material to the patient and health care worker during reuse of the transducer.

**WARNING**

- In the USA, federal law restricts device to sale or use by or on order of a physician.
- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's manual.
- Dispose of single-use components as infectious waste.
- See your system's manual for reprocessing transducer between use.



#### COVERING THE TRANSDUCER AND BRACKET

**WARNING**

- Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.
- Disposable components are packaged sterile and are single-use only. Do not use if integrity of packaging is violated or if expiration date has passed.
- For illustration purposes only, transducer may be shown without a transducer cover. Always place a cover over transducer to protect patients and users from cross-contamination.

- Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
- Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
- Secure cover with enclosed bands.
- Inspect cover to ensure there are no holes or tears.

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COPYRIGHT © 2004 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL INSTRUMENTS CO., INC. NEOGUARD IS A TRADEMARK OF CIVCO MEDICAL INSTRUMENTS CO., INC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MPS GmbH