
Laparoscopic Transducer Cover with Latex-Free Tip

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



紫外線に感光します
Sensível à luz UV
Sensible a rayos ultravioleta
Känslig mot UV-ljus

Laparoscopic Transducer Cover with Latex-Free Tip English

INTENDED USE

- *The cover allows use of the transducer in scanning and needle guided procedures for body surface, endocavity and intra-operative diagnostic ultrasound, while helping to prevent transfer of microorganisms, body fluids and particulate material to the patient and health care worker during reuse of the transducer.*

WARNING

- *In the USA, federal law restricts device to sale or use by or on order of a physician.*
- *Before use you should be trained or under the supervision of someone who is trained in ultrasonography. For instructions on preparation , safe operation and care of your transducer, see your system's operator manual.*
- *Sterile products have been sterilized by 100% ethylene oxide gas. Covers cannot be resterilized. Use only water-soluble agents.*

COVERING THE TRANSDUCER

- Inspect cover to ensure there are no holes or tears.
- Lubricate probe shaft with Endo-Lube™ to aid in application and removal of cover. Place appropriate amount of saline solution inside cover. Poor imaging may result if no coupling agent is used.

USING TRANSDUCER WITH TROCAR

- Cover can be used on 10mm transducers only. To accommodate use of transducer cover, a trocar 2mm larger than transducer size is recommended for use.
- Lubricate exterior of latex cover and trocar seal, as indicated in Endo-Lube instructions.
- To avoid tearing cover, transducer must be inserted or removed in non-angled position into trocar.
- Maintain taut fit of cover on transducer.

DISPOSAL

- Dispose of single-use components after use according to medical regulations for biohazardous waste.

CLEANING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- See your system's operator manual for instructions on cleaning and disinfecting your transducer.

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COPYRIGHT © 2005 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL INSTRUMENTS CO., INC. ULTRA-PRO II IS A TRADEMARK OF CIVCO MEDICAL INSTRUMENTS CO., INC. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA.
MPS GmbH