


Reusable Endocavity Needle Guide

*For use with Siemens 6.5EV13, 6.5EV13S, 6.5EV13+
and EV9-4 (on G20™, G50™ and G60 S™ systems)
transducers*

Reference Guide

INTENDED USE

- The reusable endocavity needle guide provides physicians with a tool for guiding a needle instrument into the ovaries or uterine wall with the use of diagnostic ultrasound transducers.

WARNING

- *In the USA, federal law restricts device to sale or use by or on order of a physician.*
- *Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's manual.*

I. COVERING THE TRANSDUCER

WARNING

- *Transducer covers containing natural rubber latex may cause allergic reactions. In the USA, refer to FDA Medical Alert MDA91-1.*
- *Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.*
- *Disposable components are packaged sterile and are single-use only. Do not use if integrity of packaging is violated or if expiration date has passed.*
- *For illustration purposes only, transducer may be shown without a transducer cover. Always place a cover over transducer to protect patients and users from cross-contamination.*
- *Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.*

1. Place an appropriate amount of water-based gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
2. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
3. Secure cover with enclosed bands.
4. Inspect cover to ensure there are no holes or tears.

II. ATTACHING NEEDLE GUIDE TO TRANSDUCER

WARNING

- *Needle guide is packaged non-sterile and is reusable. To avoid possible patient contamination, ensure bracket is properly cleaned, sterilized or disinfected before each use.*
- *Do not use needle guide if it appears damaged.*
- *Use care when inserting needle to avoid bending.*

1. With clamp open, place hook around middle of transducer and slide forward until properly secure on end of transducer.
2. Position end of needle guide over colored ring, close around transducer shaft and turn clamp screw to tighten. Check for secure fit.

NEEDLE PATH VERIFICATION

- Before using component for first time, perform needle path verification to verify system and biopsy guide relationships, as described in system manual.
- Use appropriate needle length to reach target area.

III. DISPOSAL

- Dispose of single-use components as infectious waste. Clean and sterilize reusable components after each use.

IV. CLEANING, DISINFECTING AND STERILIZING NEEDLE GUIDE

WARNING

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*
- *See your system’s manual for reprocessing transducer between use.*
- *The recommended protocol for reprocessing the reusable needle guide has been evaluated for effectiveness and compatibility. If choosing to use a reprocessing protocol other than those listed, a validation must first be completed to ensure effectiveness and compatibility.*

CLEANING THE NEEDLE GUIDE

1. After each use, remove needle guide from transducer, discarding single-use, disposable components. Remove visible contaminants from surface using a small, soft instrument brush. Keep parts from drying out until complete cleaning can be accomplished.
2. Soak needle guide for minimum of five minutes in neutral pH, low foaming enzymatic detergent, such as ENZOL® Enzymatic Detergent (Johnson & Johnson). While immersed, use instrument brush to remove trapped contaminants from surfaces and features. If visible contaminants cannot be easily removed, repeat soaking procedure for an additional five minutes. Remove needle guide from cleaning solution and remove any remaining residue with dry wipe. Follow cleaning solution manufacturer’s directions for use and recommendations for concentration.

DISINFECTING AND STERILIZING NEEDLE GUIDE

- High-level disinfect needle guide using CIDEX OPA® ortho-Phthalaldehyde (Johnson & Johnson) or equivalent .55% ortho-phthalaldehyde-based solution. In the U.S.A., use only disinfecting/sterilizing agent cleared by FDA and approved by EPA, such as CIDEX® Activated Dialdehyde Solution. Follow manufacturer’s instructions and recommendations for concentration, time of contact and post-process procedure.
- High-level disinfect or sterilize needle guide using CIDEX® Activated Dialdehyde (Johnson & Johnson) or equivalent 2% glutaraldehyde-based solution, CIDEX Plus® (Johnson & Johnson) or equivalent 3.4% glutaraldehyde-based solution, or a hydrogen peroxide-based solution. Needle guide may be sterilized by using steam (autoclave) at 132-138°C for 15 to 30 minutes.
- Inspect parts for damage such as cracks or breakage. If damage is evident, discontinue use of needle guide and contact your ultrasound system support representative.

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COPYRIGHT © 2012. ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
MPS GmbH