# Reusable Endocavity Needle Guide

---

REFERENCE GUIDE

CE  
0120

**I. Covering the Transducer**

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

**INTENDED USE**

- This device provides physicians with a tool for performing needle/instrument guided procedures with the use of the diagnostic ultrasound endocavity transducer.

**WARNING**

- *In the USA, federal law restricts device to sale or use by or on order of a physician.*
- *Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system’s manual.*
- *For illustration purposes only, transducer may be shown without a transducer cover. Always place a cover over transducer to protect patients and users from cross-contamination.*
- *Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.*
- *Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.*

**I. Covering the Transducer**

1. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
2. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
3. Secure cover with enclosed bands.
4. Inspect cover to ensure there are no holes or tears.

**II. Attaching Needle Guide to Transducer**

**WARNING**

- *Needle guide is packaged non-sterile and is reusable. To avoid possible patient contamination, ensure needle guide is properly cleaned, sterilized or disinfected before each use.*
  - *Do not use needle guide if it appears damaged.*
1. With clamp open, align guide with locating features on transducer.
  2. Close needle guide around transducer and check for secure fit.

**NEEDLE PATH VERIFICATION**

3. Before using component for first time, perform needle path verification to verify system and biopsy guide relationships, as described in system manual.
4. Use appropriate needle length to reach target area.

**III. Disposal**

**WARNING**

- *Dispose of single-use components as infectious waste.*

**IV. Cleaning, Disinfecting and Sterilizing Needle Guide**

**WARNING**

- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*
- *See your system’s manual for reprocessing transducer between use.*
- *The recommended protocol for reprocessing the reusable needle guide has been evaluated for effectiveness and compatibility. If choosing to use a reprocessing protocol other than those listed, a validation must first be completed to ensure effectiveness and compatibility.*

---

## *CLEANING THE NEEDLE GUIDE*

---

1. After each use, remove needle guide from transducer, discarding single-use, disposable components. Remove visible contaminants from needle guide surface using a small, soft instrument brush. Keep needle guide from drying out until complete cleaning can be accomplished.
2. Soak needle guide for minimum of five minutes in neutral pH, low foaming enzymatic detergent, such as ENZOL<sup>®</sup> Enzymatic Detergent (Johnson & Johnson). While immersed, use instrument brush to remove trapped contaminants from surfaces and features. If visible contaminants cannot be easily removed, repeat soaking procedure for an additional five minutes. Remove needle guide from cleaning solution and remove any remaining residue with dry wipe. Follow cleaning solution manufacturer's directions for use and recommendations for concentration.

---

## *DISINFECTING AND STERILIZING NEEDLE GUIDE*

---

- High-level disinfect needle guide using CIDEX OPA<sup>®</sup> ortho-Phthalaldehyde Solution (Johnson & Johnson) or equivalent .55% ortho-phthalaldehyde-based solution. Follow manufacturer's instructions and recommendations for concentration, time of contact and post-process procedure.
- High-level disinfect or sterilize needle guide using CIDEX<sup>®</sup> Activated Dialdehyde Solution (Johnson & Johnson) or equivalent 2% glutaraldehyde-based solution, CIDEX Plus<sup>®</sup> (Johnson & Johnson) or equivalent 3.4% glutaraldehyde-based solution, or a hydrogen peroxide-based solution. Needle guide may be sterilized by using steam (autoclave) at 132-138°C for 15 to 30 minutes.
- Inspect needle guide for damage such as cracks or breakage. If damage is evident, discontinue use of needle guide and contact your ultrasound system support representative.

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COPYRIGHT © 2016 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
MPS GmbH