DISPOSABLE ENDOCAVITY NEEDLE GUIDE



For use with Exact Imaging, Mindray, and SuperSonic Imagine transducers

REFERENCE GUIDE



CE
1639

Exact Imaging EV9C



16-18 GA

Mindray 65EB10EA, V10-4s, V11-3s, V11-3Ws

16-18 GA

Supersonic Imagine E12-3, SE12-3, SEV12-3

16-18 GA

NEEDLE GUIDE

INDICATIONS FOR USE

Diagnostic ultrasound needle / instrument guided procedures - transvaginal and transrectal: tissue biopsy, fluid aspiration, catheter placement, and treatment.

INTENDED USE

This device provides a mechanical means for performing needle / instrument guided procedures with the use of the diagnostic ultrasound endocavity transducer. Needle guide is attached over endocavity transducer / probe / scanhead instruments. This device provides a fixed path for the needle or instrument that when coupled by the ultrasound system corresponds to on-screen imaging guidelines for visualizing guided instrument placement procedures. CIVCO Endocavity Ultrasound Needle / Biopsy / Guides are furnished sterile; single use patient / procedure, disposable. The single use, disposable feature helps to prevent transfer of microorganisms, body fluids, and particle material to the patient and healthcare worker during reuse of the transducer.

COVER

INDICATIONS FOR USE / INTENDED USE

Protective cover or sheath placed over diagnostic ultrasound transducer / probe / scanhead instruments. The cover allows use of the transducer in scanning and needle guided procedures for body surface, endocavity, and intra-operative diagnostic ultrasound, while helping to prevent transfer of microorganisms, body fluids, and particulate material to the patient and healthcare worker during reuse of the transducer (both sterile and non-sterile covers). The cover also provides a means for maintenance of a sterile field (sterile covers only). CIVCO Latex Transducer Covers are furnished sterile & non-sterile; single use patient / procedure, disposable.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's user guide.
- Transducer covers containing natural rubber latex may cause allergic reactions. In the USA, refer to FDA Medical Alert MDA91-1.
- Use of a needle guide on an endocavity transducer may increase the risk of vaginal or rectal abrasion.
- For illustration purposes only, transducer may be shown without a transducer cover. Always place a cover over transducer to protect patients and user from cross-contamination.
- Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.
- Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.
- Disposable components are packaged sterile and are single-use only. Do not use if integrity of packaging is violated or if expiration date has passed.
- Do not use if needle guide is damaged or does not fit properly.

NEEDLE PATH VERIFICATION


- Before using component for first time, perform needle path verification to verify system and biopsy guide relationships, as described in system user guide.
- See transducer user guide for proper procedure.
- Use care when inserting needle to avoid bending.

NOTE: Use appropriate needle length to reach target area.


USING NEEDLE GUIDANCE

1. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
2. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
3. Secure cover with enclosed bands.


4. Inspect cover to ensure there are no holes or tears.


5. Using proper sterile technique, align guide with locating features of transducer.

6. Press down to secure to transducer.

7. If desired, apply secondary cover by following the steps above.

DISPOSAL

- *Dispose of single-use components as infectious waste.*
- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*
- *See your system's user guide for reprocessing transducer between use.*

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COPYRIGHT © 2020 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
MPS GmbH
CIVCO Medical Instruments Co., Inc.