Disposable Transrectal Needle Guide


For use with GE Medical Systems
7.5 Endorectal transducer

REFERENCE GUIDE


CE
0120


GE Medical Systems 7.5 Endorectal (Ausonics)

1

2

16-18GA

3

INTENDED USE

- This device provides physicians with a tool for performing needle/instrument guided procedures with the use of the diagnostic ultrasound endocavity transducer.

WARNING

- *In the USA, federal law restricts device to sale or use by or on order of a physician.*
- *Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's manual.*
- *Transducer covers containing natural rubber latex may cause allergic reactions. In the USA, refer to FDA Medical Alert MDA91-I.*
- *Disposable components are packaged sterile and are single-use only. Do not use if integrity of packaging is violated or if expiration date has passed.*
- *For illustration purposes only, transducer may be shown without a transducer cover. Always place a cover over transducer to protect patients and users from cross-contamination.*
- *Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.*
- *Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.*

COVERING THE TRANSDUCER

1. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.
2. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
3. Secure cover with enclosed bands.
4. Inspect cover to ensure there are no holes or tears.

ATTACHING NEEDLE GUIDE TO TRANSDUCER

1. Using proper sterile technique, align guide with locating features of transducer.
2. Align guide with locating features on transducer.
3. Press down to secure to transducer, locating the slide bar slot on transducer.
4. To scan, slide needle guide backward until it stops.
5. To biopsy, slide needle guide forward until it stops.

NEEDLE PATH VERIFICATION

WARNING

- *Before using component for first time, perform needle path verification to verify system and biopsy guide relationships.*
- *Use appropriate needle length to reach target area.*
- *See transducer manual for proper procedure.*
- *Use care when inserting needle to avoid bending.*

DISPOSAL

WARNING

- *Dispose of single-use components as infectious waste.*
- *Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.*
- *See your system's manual for reprocessing transducer between use.*

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COPYRIGHT © 2012 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
MPS GmbH