ACCUSITE™

NEEDLE GUIDANCE SYSTEM

REFERENCE GUIDE

For use with
Canon Medical Systems
PLI-2004BX (i24LX8) and PLU-704SDT
transducers

CE
0120

BRACKET / NEEDLE GUIDE

INDICATIONS FOR USE

The guidance system is intended for directing instruments such as catheters, electrodes and needles into a targeted anatomical location of a patient relative to the imaging instrument for percutaneous procedures. Prior to use healthcare workers should be trained in ultrasonography.

INTENDED USE

The bracket and disposable guide provide physicians with a tool for performing needle-guided (or catheter) procedures with the use of diagnostic ultrasound transducers.

COVER

INDICATIONS FOR USE / INTENDED USE

Protective cover or sheath placed over diagnostic ultrasound transducer / probe / scanhead instruments. The cover allows use of the transducer in scanning and needle guided procedures for body surface, endocavity, and intra-operative diagnostic ultrasound, while helping to prevent transfer of microorganisms, body fluids, and particulate material to the patient and healthcare worker during reuse of the transducer (both sterile and non-sterile covers). The cover also provides a means for maintenance of a sterile field (sterile covers only). CIVCO Poly Ultrasound Transducer Covers are furnished sterile & non-sterile; single use patient / procedure, disposable.

CAUTION

Federal (United States) law restricts this device to sale by or on the order of a physician.

WARNING

- Before use, you should be trained in ultrasonography. For instructions on the use of your transducer, see your system's user guide.
- Disposable components are packaged sterile and are single-use only. Do not use if integrity of packaging is violated or if expiration date has passed.
- Do not reuse, reprocess or resterilize single-use device. Reuse, reprocessing or resterilization may create a risk of contamination of the device, cause patient infection or cross-infection.
- CIVCO's needle guidance systems are designed and validated as an integrated two-part system consisting of a custom reusable bracket and a disposable needle guide. Use of non-CIVCO guides on a CIVCO bracket may result in serious injury to the patient and/or user. To ensure safe and effective use in compliance with system software guidelines, only CIVCO needle guides are recommended.
- Bracket is packaged non-sterile and is reusable. To avoid possible patient contamination, ensure bracket is properly cleaned, sterilized or disinfected before each use.
- Do not use if bracket is damaged or does not fit properly.
- Do not use if needle guide is damaged or does not fit properly.
- Use only water-soluble agents or gels. Petroleum or mineral oil based materials may harm cover.
- For illustration purposes only, transducer and bracket may be shown without a transducer cover. Always place a cover over transducer and bracket to protect patients and users from cross-contamination.


NEEDLE PATH VERIFICATION

- Before using component for first time, perform needle path verification to verify system and biopsy guide relationships, as described in system user guide.

NOTE: Use appropriate needle length to reach target area.


USING NEEDLE GUIDANCE SYSTEM

1. Attach bracket to transducer by aligning locating features. Ensure bracket is firmly attached.


2. Place an appropriate amount of gel inside cover and/or on transducer face. Poor imaging may result if no gel is used.

3. Insert transducer into cover making sure to use proper sterile technique. Pull cover tightly over transducer face to remove wrinkles and air bubbles, taking care to avoid puncturing cover.
4. Secure cover with enclosed bands.
5. Inspect cover to ensure there are no holes or tears.


6. Select the appropriate needle guide to achieve target depth from skin line.
7. Using proper sterile technique, hook needle guide onto bracket and snap onto attachment area, taking care to avoid puncturing cover.
8. Ensure the door on needle guide is fully closed.


QUICK RELEASE FUNCTION

1. Depress tab on needle guide toward bracket to activate needle guide quick release function.

REMOVING NEEDLE GUIDE

1. Press tab away from bracket to remove needle guide.

CLEANING, DISINFECTING AND STERILIZING BRACKET

WARNING

- Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.
- Dispose of single-use components as infectious waste. Clean and sterilize reusable components after each use.
- See your system's user guide for reprocessing transducer between use.
- This procedure for reprocessing has been validated for effectiveness and compatibility. Device could be damaged or cross-contaminated due to improper reprocessing.
- Do not gas sterilize or autoclave bracket.
- Inspect bracket for damage such as cracks or breakage. If damage is evident, discontinue use of bracket and contact your ultrasound system support representative.

CLEANING THE BRACKET

1. After each use, remove bracket and needle guide from transducer, discarding single-use, disposable components. Remove visible contaminants from bracket surface using a small, soft instrument brush. Keep bracket from drying out until complete cleaning can be accomplished.
2. Soak bracket for minimum of five minutes in neutral pH, low foaming enzymatic detergent, such as ENZOL® Enzymatic Detergent (Johnson & Johnson). While immersed, use instrument brush to remove trapped contaminants from surfaces and features. If visible contaminants cannot be easily removed, repeat soaking procedure for an additional five minutes. Remove bracket from cleaning solution and remove any remaining residue with dry wipe. Follow cleaning solution manufacturer's directions for use and recommendations for concentration.

DISINFECTING AND STERILIZING BRACKET

- High-level disinfect bracket using CIDEX OPA® ortho-Phthalaldehyde Solution (Johnson & Johnson) or equivalent .55% ortho-phthalaldehyde-based solution. Follow manufacturer's instructions and recommendations for concentration, time of contact and post-process procedure.
- High-level disinfect or sterilize bracket using CIDEX® Activated Dialdehyde Solution (Johnson & Johnson) or equivalent 2% glutaraldehyde-based solution, CIDEX Plus® (Johnson & Johnson) or equivalent 3.4% glutaraldehyde-based solution, or a hydrogen peroxide-based solution. Follow manufacturer's instructions and recommendations for concentration, time of contact and post-process procedure.

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COPYRIGHT © 2018 ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ACCUSITE IS A TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL OTHER TRADEMARKS ARE PROPERTY OF THEIR RESPECTIVE OWNERS. PRINTED IN USA. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW.
MPS GmbH
CIVCO Medical Instruments Co., Inc.