

ASTRA TEE® QUICK START INSTRUCTIONS

	At rest, the display reads:	<div style="border: 1px solid black; border-radius: 15px; padding: 10px; width: fit-content; margin: 0 auto;"> <p>Ready for Use Press Enter to Start</p> </div>
Step 1:	Press ENTER. Display will read:	<div style="border: 1px solid black; border-radius: 15px; padding: 10px; width: fit-content; margin: 0 auto;"> <p>Scan Probe</p> </div>
Step 2:	Scan barcode on probe electrical connector and place into the wall-mounted Connector Holder. Display will read:	<div style="border: 1px solid black; border-radius: 15px; padding: 10px; width: fit-content; margin: 0 auto;"> <p>Insert Probe</p> </div>
Step 3:	Open door and insert pre-cleaned probe into disinfection chamber, placing probe handle into Probe Handle Holder. Loop probe cable over cable strain relief. Display will read:	<div style="border: 1px solid black; border-radius: 15px; padding: 10px; width: fit-content; margin: 0 auto;"> <p>Scan Exam ID ENTER = skip</p> </div>
Step 4:	Entering Exam ID is optional, but you must scan ID or press ENTER to continue. Display will read:	<div style="border: 1px solid black; border-radius: 15px; padding: 10px; width: fit-content; margin: 0 auto;"> <p>Close Door or Scan Second Probe</p> </div>
Step 5:	Install second probe (if required) following same instructions or close door. Barcode scanner will remain illuminated until door is closed or a second probe is scanned. Display will read:	<div style="border: 1px solid black; border-radius: 15px; padding: 10px; width: fit-content; margin: 0 auto;"> <p>Scan Operator ID ENTER = Skip</p> </div>

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

<p>Step 6:</p>	<p>Scanning Operator ID or pressing ENTER to Skip will continue to the next step of scanning the test strip bottle. Display will read:</p>	<p style="text-align: center;">Scan TS Bottle MENU = Cancel</p>
<p>Step 7:</p>	<p>As soon as the test strip bottle is scanned, HLD is then pumped up into the disinfection chamber. Display will read:</p>	<p style="text-align: center;">Filling chamber with HLD, please wait....</p>
<p>Step 8:</p>	<p>Once the fill level has been reached, the MRC process may begin. Display will read:</p>	<p style="text-align: center;">MRC Countdown Timer Skip < -- > Start</p>
<p>Step 9:</p>	<p>Utilizing the MRC Countdown Timer is optional, but you must press either Skip or Start to continue. If Skip is selected continue to Step 10. If you plan to use the timer, open disinfection chamber door. Before selecting Start, ensure test strip is in hand, ready to use. Select Start and dip an MRC test strip into HLD, holding it submerged for dip time. Display will read:</p>	<p style="text-align: center;">TS DIP READ XXs XXs</p>
<p>Step 10:</p>	<p>Once the countdown timer is complete (or skipped), the results of the MRC test may be read and documented. Display will read:</p>	<p style="text-align: center;">Test MRC PASS < > FAIL</p>
<p>Step 11:</p>	<p>Compare test strip to color chart on manufacturer's bottle. Display will read:</p> <p>Press PASS or FAIL using the corresponding arrows.</p> <p>NOTE: Pressing FAIL requires a new bottle of HLD.</p>	<p style="text-align: center;">Test MRC PASS < > FAIL</p>
<p>Step 12:</p>	<p>If PASS is pressed, display will read:</p> <p>ONE PROBE: If desired, conduct leak testing according to manufacturer's instructions. Press PASS or FAIL using corresponding arrows. Data will be logged.</p> <p>TWO PROBES: If desired, selecting which probe to test first using the corresponding arrows. Conduct leak testing according to manufacturer's instructions. Press PASS or FAIL using the corresponding arrows. Repeat process with other probe. Data will be logged.</p>	<p style="text-align: center;">Test Probe Leakage PASS <- -> FAIL</p> <p style="text-align: center;">Probe to test first? Left <- -> Right</p>

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Step 13:	Display will read:	
Step 14:	On successful completion of disinfection and rinse cycles, display will read:	
Step 15:	The probe(s) are disinfected and should be removed. After removal, display will read:	

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