



ULTRASOUND IMAGING ACCESSORIES

REGIONAL ANESTHESIA



VOLUME 1
WWW.CIVCO.COM





REGIONAL ANESTHESIA

Welcome to the 2010 CIVCO Medical Solutions Regional Anesthesia catalog, part of our new Ultrasound Imaging Accessories catalog collection.

We are pleased you have chosen CIVCO's high-quality products and look forward to serving you. Our catalog will provide you with quick and easy access to the latest in regional anesthesia products for use with leading ultrasound equipment.

This year we divided our Ultrasound Imaging Accessories catalog into a collection of individualized specialty catalogs in order to better serve our customers. By offering products that meet physicians' specific procedural needs, we hope to make the ordering process quicker and easier.

As a valued customer, CIVCO is committed to providing you with innovative and clinically beneficial products and services. Special care is taken for each and every product going out our door.

For over two decades, CIVCO has been providing imaging accessories to facilities around the world. Our company goal remains the same today, "Total Customer Satisfaction." If you are not completely satisfied with any CIVCO product, please contact one of our customer service specialists and they will be happy to assist you.

Thank you for your business. We look forward to serving you again in the near future.



DIRECTORY

<u>General Purpose Transducer Covers</u>	<u>6</u>
<u>General Purpose Needle Guides</u>	<u>7</u>
<u>OEM Partner Accessory Guide</u>	<u>9</u>
<u>General Purpose Accessories</u>	<u>10</u>
<u>Specialty Covers & Markers</u>	<u>11</u>
<u>Transmission Gel</u>	<u>12</u>
<u>Cleaners & Disinfectants</u>	<u>12</u>
<u>Ordering Information</u>	<u>15</u>
<u>Index By Part Number</u>	<u>16</u>



One Company. Many Solutions.

Featured Solutions

Extended Length CIV-FLEX™ Covers

Providing patient and staff safety.

Recommended for over two decades as the latex-free cover of choice. The improved telescopic-fold allows for easy application of gel, as well as easy application of the cover onto the transducer.

FEATURES & BENEFITS

- Extended length covers offer additional sterile protection
- Available in extended lengths of 36 and 58 inches
- Select covers offer a three-dimensional "box end"

See page 6 for more details or call 319.656.4447 for pricing



AccuSITE™ Needle Guide

Safe. Efficient. Guiding Patient Care.

The AccuSITE needle guide is designed specifically for guidance in central line placement procedures. Designed for minimizing hospital risks and increasing patient well being.

FEATURES & BENEFITS

- First time success rates for accuracy
- Features large funnel for easy needle insertion
- Disposable, single-use guide increases infection prevention
- Available in 18, 20 and 21/22 gauge
- Depths include: .5, 1, 1.5, 2, 2.5, 3 and 3.5cm
- Features colored tab for quick-release function

See page 8 for more details or call 319.656.4447 for pricing



Infiniti™ Needle Guide

Providing accurate in-plane guidance.

The Infiniti guide offers multiple angle capabilities for shallow and deep access. The open channel design of the guide allows instruments to be positioned in the scan plane, providing a safe, efficient means for ultrasound-guided puncture procedures.

FEATURES & BENEFITS

- Single-use design increases infection prevention
- Design allows users an efficient needle release for viewing anatomy in the long axis
- Accepts 14 and 18 or 20, 21/22 and 25 gauge instruments

See page 9 for more details or call 319.656.4447 for pricing



SonoSite L25 Series Needle Guide

Providing clinical confidence for line placement and nerve blocks.

The SonoSite L25 Series needle guide provides a fixed path for visualizing ultrasound-guided vascular access and nerve block procedures. The guide features a mechanism for quick-release detachment of the needle from the guide while maintaining secure needle placement once the targeted depth has been reached.

FEATURES & BENEFITS

- Single-use design reduces risks of cross-contamination
- Snaps onto bracket with single-click, no additional assembly is required
- Colored side tab allows for simple quick-release function
- Available in 18, 21 and 22 gauge
- Depths include: 1, 1.5, 2, 2.5, 3.5 and 4.5cm



NOW AVAILABLE DIRECT FROM CIVCO! See page 8 for more details or call 319.656.4447 for pricing



Latex-Free Telescopically-Folded Extended-Length CIV-Flex™ Covers

- Extended lengths of 36" and 58" offer additional sterile protection
- Sterile procedure kits include: transducer cover, sterile gel packet and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
610-637	Sterile 8.9 x 91.5cm (3.5" x 36") CIV-Flex cover	24
610-1000	Sterile 10.2 x 147cm (4" x 58") CIV-Flex cover	24
610-1064	Sterile 10.2 x 147cm (4" x 58") CIV-Flex cover (3D) (no gel)	24
610-542	Sterile 14 x 91.5cm (5.5" x 36") CIV-Flex cover (3D)	24
610-575	Sterile 14 x 147cm (5.5" x 58") CIV-Flex cover (3D)	24
610-1010	Sterile 17.8 x 147cm (7" x 58") CIV-Flex cover (3D)	24



Latex-Free Flat-Folded CIV-Flex™ Covers

- Soft, pliable latex-free material
- Offers sterile protection when performing routine puncture and drainage procedures
- Sterile procedure kits include: transducer cover, sterile gel packet and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
610-323	Sterile 6 x 91.5cm (2.5" x 36") CIV-Flex cover	24
610-001	Sterile 8.9 x 61cm (3.5" x 24") CIV-Flex cover	24
610-002	Sterile 14 x 61cm (5.5" x 24") CIV-Flex cover (3D)	24
610-004 (04510947)	Sterile 17.8 x 61cm (7" x 24") CIV-Flex cover	24
610-005	Sterile 22.9 x 61cm (9" x 24") CIV-Flex cover	24



Latex-Free iVAS™ Vascular Access Transducer Cover

- Recommended for vascular access procedures including PICC, CVC and femoral artery placement
- Expands sterile field up to 58" (147cm)
- Telescopically-folded and easily opened
- Covers packaged with sterile gel packet and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
610-989	Sterile 15.2 x 147cm (6" x 58") telescopically-folded polyethylene cover	24



Latex-Free Polyethylene Covers

- Economically-priced
- Safe means for reducing cross-contamination
- Select covers offer 3D "box end"
- Covers packaged with sterile gel packet and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
610-1050	Sterile 15.2 x 25.4cm (6" x 10") polyethylene cover (no gel)	24
610-013	Sterile 5 x 61cm (2" x 24") polyethylene cover	24
610-015	Sterile 10.2 x 61cm (4" x 24") polyethylene cover (3D)	24
610-208	Sterile 15.2 x 83.8cm (6" x 33") polyethylene cover	24
610-020	Sterile 17.8 x 61cm (7" x 24") polyethylene cover	24
610-022	Sterile 22.9 x 61cm (9" x 24") polyethylene cover	24



Latex-Free NeoGuard™ Covers

- Rolled, latex-free material specially designed for ultrasound use
- Recommended for patients and professionals with latex protein sensitivities
- Available for small footprint transducers
- Each cover includes colored elastic bands
- Sterile gel packets available separately

REF. #	PRODUCT DESCRIPTION	QTY.
610-844	Sterile 4 x 30cm (1.6" x 11.8") NeoGuard cover	24
610-838	Non-sterile 4 x 30cm (1.6" x 11.8") NeoGuard cover	24
610-1038	Sterile 2.6 x 30cm (1" x 11.8") NeoGuard cover	50



Latex Covers

- Rolled for easy application to a variety of transducers
- Extended lengths available
- Offers sterile coverage during general puncture and drainage procedures
- Covers packaged with sterile gel packet and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
610-042	Sterile 4 x 30cm (1.6" x 11.8") latex cover	24
610-043	Sterile 5 x 30cm (2" x 11.8") latex cover	24
610-044	Sterile 8 x 45cm (3.2" x 18") latex cover	24
610-046	Sterile 10 x 45cm (3.9" x 18") latex cover	24



NEW!

Ultra-Pro e™ Disposable Variable Angle Needle Guides

- Open channel offers infinite angle capabilities
- Ensures accurate in-plane positioning during puncture procedures
- Accepts 14GA and smaller instruments
- Designed for use with Ultra-Pro™ or Ultra-Pro II™ style transducer brackets
- U.S. Patent No. 6,379,307

CHECK WEBSITE FOR AVAILABILITY

REF. #	PRODUCT DESCRIPTION	QTY.
610-1088	Sterile Ultra-Pro e guide, 14GA-green	24
610-1089	Sterile Ultra-Pro e guide, 14GA-green, with 14 x 91.5cm (5.5" x 36") telescopically-folded CIV-Flex cover	24
610-1090	Sterile Ultra-Pro e guide, 14GA-green, with 14 x 147cm (5.5" x 58") telescopically-folded CIV-Flex cover	24
610-1082	Sterile Ultra-Pro e guide, 16GA-white	24
610-1083	Sterile Ultra-Pro e guide, 16GA-white, with 14 x 91.5cm (5.5" x 36") telescopically-folded CIV-Flex cover	24
610-1084	Sterile Ultra-Pro e guide, 16GA-white, with 14 x 147cm (5.5" x 58") telescopically-folded CIV-Flex cover	24

Visit WWW.CIVCO.COM for photos, pricing and system spec sheets.



AccuSITE™ Disposable Transverse Needle Guides

- Increases user productivity during central line placement procedures
- Intuitive, easy-to-use design
- Single-click assembly
- Short axis guide provides a safe, reliable and efficient tool for vessel cannulation
- Large, colored quick-release tab
- Accepts 18, 20 and 21/22GA instruments
- Depths include: .5, 1, 1.5, 2, 2.5, 3 and 3.5cm
- U.S. Patents Pending

REF. #	PRODUCT DESCRIPTION	QTY.
610-1017	Sterile AccuSITE needle guide with 14 x 147cm (5.5" x 58") CIV-Flex cover (3D)-18 gauge only. (Pink)	24
610-1030	Sterile AccuSITE needle guide with 14 x 147cm (5.5" x 58") CIV-Flex cover (3D)-18 gauge only. (Pink)	Demo Kit
610-1018	Sterile AccuSITE needle guide with 14 x 147cm (5.5" x 58") CIV-Flex cover (3D)-20 gauge only. (Yellow)	24
610-1031	Sterile AccuSITE needle guide with 14 x 147cm (5.5" x 58") CIV-Flex cover (3D)-20 gauge only. (Yellow)	Demo Kit
610-1019	Sterile AccuSITE needle guide with 14 x 147cm (5.5" x 58") CIV-Flex cover (3D)-21/22 gauge only. (Green)	24
610-1032	Sterile AccuSITE needle guide with 14 x 147cm (5.5" x 58") CIV-Flex cover (3D)-21/22 gauge only. (Green)	Demo Kit



Now Available from CIVCO!

SonoSite L25 series Needle Guides

For use with SonoSite L25 series transducers

- L25 needle guide allows for accurate short axis placement
- Quick-release mechanism for easy detachment of needle from guide
- Accepts 18, 21 and 22GA instruments
- Depths include: 1, 1.5, 2, 2.5, 3.5 and 4.5cm
- U.S. Patent No. 7,588,541

REF. #	PRODUCT DESCRIPTION	QTY.
P09072 (P06014)	Sterile needle guides (1.5cm, 2.5cm, 3.5cm, 4.5cm depths) with 10.2 x 147cm (4" x 58") telescopically-folded CIV-Flex cover, gel packet and bands. 18 gauge only. (Blue)	24
P09073 (P06015)	Sterile needle guides (1cm, 1.5cm, 2cm depths) with 10.2 x 147cm (4" x 58") telescopically-folded CIV-Flex cover, gel packet and bands. 21 gauge only. (Teal)	24
P09074 (P06016)	Sterile needle guides (1cm, 1.5cm, 2cm depths) with 10.2 x 147cm (4" x 58") telescopically-folded CIV-Flex cover, gel packet and bands. 22 gauge only. (Black)	24

L25 Series Needle Guide is Now Available Directly from CIVCO!



NEW! Smaller Gauge Guides

Infiniti™ Disposable Needle Guides

- Offers multiple angle capabilities for shallow and deep access
- Open channel design of needle guide allows instruments to be positioned in the scan plane
- Accepts 14 and 18GA or 20, 21/22 and 25GA instruments
- Replacement kits include: needle guide, latex-free transducer cover, gel packet and colored elastic bands
- U.S. Patent No. 6,379,307

CHECK WEBSITE FOR AVAILABILITY

REF. #	PRODUCT DESCRIPTION	QTY.
610-1073	Sterile Infiniti needle guide-20, 21/22, 25GA	24
610-1074	Sterile Infiniti needle guide-20, 21/22, 25GA with 14 x 91.5cm (5.5" x 36") CIV-Flex cover (3D)	24
610-1075	Sterile Infiniti needle guide-20, 21/22, 25GA with 14 x 147cm (5.5" x 58") CIV-Flex cover (3D)	24
610-1094	Sterile Infiniti needle guide-14, 18GA	24
610-1095	Sterile Infiniti needle guide-14, 18GA with 14 x 91.5cm (5.5" x 36") CIV-Flex cover (3D)	24
610-889	Sterile Infiniti needle guide with 14 x 147cm (5.5" x 58") telescopically-folded CIV-Flex cover (3D)	24

Ultrasound Guided Regional Anesthesia Needle Guidance Systems

SONOSITE Transducers	Biopsy Depth	General Purpose Starter Kits
C60 series	N/A	674-036
L25 series	N/A	674-030
L25 series	N/A	674-037 (P12920)-18GA 674-038 (P12922)-21GA 674-039 (P12924)-22GA
L38 series	N/A	674-031

GE HEALTHCARE Transducers	Biopsy Depth	General Purpose Starter Kits
12L-RS	.5, 1, 1.5, 2, 2.5, 3, 3.5cm	742-359
12L-RS	N/A	742-397
12L-SC	.5, 1, 1.5, 2, 2.5, 3, 3.5cm	742-359
12L-SC	N/A	742-397

PHILIPS Transducers	Biopsy Depth	General Purpose Starter Kits
L12-5	N/A	613-209
L17-5	N/A	613-247

ARROW-TELEFLEX Transducers	Biopsy Depth	General Purpose Starter Kits
InView	N/A	INV-90001-BK

TERASON Transducers	Biopsy Depth	General Purpose Starter Kits
12L5V	.5, 1, 1.5, 2, 2.5, 3, 3.5cm	683-002

ULTRASONIX Transducers	Biopsy Depth	General Purpose Starter Kits
L9-4	.5, 1, 1.5, 2, 2.5, 3, 3.5cm	684-006
L14-5	.5, 1, 1.5, 2, 2.5, 3, 3.5cm	684-006

ZONARE Transducers	Biopsy Depth	General Purpose Starter Kits
C4-1	N/A	698-013
C9-3	N/A	698-009
L8-3	.5, 1, 1.5, 2, 2.5, 3, 3.5cm	698-005
L10-5	.5, 1, 1.5, 2, 2.5, 3, 3.5cm	698-008
L14-5sp	.5, 1, 1.5, 2, 2.5, 3, 3.5cm	698-006
L14-5w	.5, 1, 1.5, 2, 2.5, 3, 3.5cm	698-012
L14-5w	N/A	698-007

ESAOTE Transducers	Biopsy Depth	General Purpose Starter Kits

Visit CIVCO's website for new Esaote releases at WWW.CIVCO.COM

Visit CIVCO's website for new product releases at WWW.CIVCO.COM OR call 319.656.4447 to speak with a Customer Service Specialist to see what solution we have for you!



Ultrasound Transducer Holder

For use with the CIVCO Assist™ Positioning Arms

- Positions transducer in a fixed location for an extended period of time
- Reduces risks of repetitive stress injuries
- Autoclavable holder is ideal for performing diagnostic and therapeutic procedures
- U.S. Patent No. 7,395,563

REF. #	PRODUCT DESCRIPTION	QTY.
810-203	Ultrasound Transducer Holder- Medium Handle	1
810-205	Non-sterile Positioning Arm EL, 40.4" (102.6cm)	1
810-200	Non-sterile Positioning Arm, 32.4" (82.3cm)	1
810-202	Multi-Purpose Holder	1
810-101	Rail Adapter, 1.1" x .38" (2.8 x 1cm)	1
810-100	Table Adapter, 1" (2.5cm)	1
610-932	Sterile 25.4 tapered to 5.1 x 121.9cm (10" tapered to 2" x 48") telescopically-folded CIV-Flex cover with Quick-Mount™ perforation	48
610-946	Sterile 25.4 tapered to 5.1 x 147cm (10" tapered to 2" x 58") telescopically-folded CIV-Flex positioning arm cover	24



Arm Supports

For use with the CIVCO Assist™ Positioning Arms

- Positions and restricts the motion of a patients arm during a procedure
- Provides anatomic support for the patients forearm or arm and may help prevent brachial plexus injuries
- U.S. Patent No. 7,395,563

REF. #	PRODUCT DESCRIPTION	QTY.
810-206	Arm Supports	2



TEE Transducer Holder

For use with Acuson V5Ms, V7M, V7B, MPT7-4, Aloka UST-5293-5, UST-52116, UST-52104, UST-52110, UST-52111, GE Healthcare 6T, 6T-RS, 9T, 9T-RS, Philips microTEE, Omni III, S7-2t, S7-3t, T6207, X7-2t, SonoSite TEE Series and ZONARE T8-3 transducers

- For use in Cardiac Anesthesia
- Attaches to most tables and surgical rails
- Reduces risk of stretching corners of patient's mouth
- Allows for stable & accurate positioning
- Ergonomic design reduces risk of injury

REF. #	PRODUCT DESCRIPTION	QTY.
811-050	TEE Holder Kit with TEE Holder (1) Rail Adapter (1) and non-sterile 4.6 tapered to 2.5 x 122cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator, gel-filled syringe and bite block (12)	1
811-051	TEE Holder Kit with TEE Holder (1), Table Adapter (1), and non-sterile 4.6 tapered to 2.5 x 122 cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator, gel-filled syringe, and bite block (12)	1
811-052	TEE Holder Kit with TEE Holder(1) Rail Adapter(1) and non-sterile 4.6 tapered to 2.5 x 122cm (1.8" tapered to 1" x 48") CIV-Flex cover on applicator, gel filled syringe and bite block(12) - Philips(Ref #8104)	1
811-053	TEE Holder Kit with TEE Holder(1), Table Adapter(1), and non-sterile 4.6 tapered to 2.5 x 122cm(1.8" tapered to 1"x48") CIV-Flex cover on applicator, gel filled syringe, and bite block(12) - Philips (Ref#8104)	1



Latex-Free Surgi-Boot™ Intraoperative Covers

For use with Biosound Esaote Omni, GE Healthcare i12L, Philips (ATL) CL10-5, CL15-7, Siemens VF13-5SP, SonoSite HST and SLA/13-6 Series

- Extends field of sterile protection against cross-contamination
- Latex-free, form-fitting tip provides snug fit to small footprint transducers
- Individually packaged with sterile gel packet and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
610-797	Sterile 15.2 x 244cm (6" x 96") poly telescopically-folded cover with attached form-fitting latex-free NeoGuard tip	12



Latex-Free Surgi-Tip™ Intraoperative Covers

- Provides snug fit to small footprint transducers
- Accordion or telescopically-folded for easy application
- Extended sterile coverage for transducer and cord
- Covers include: sterile gel packet, probe clips and colored elastic bands

REF. #	PRODUCT DESCRIPTION	QTY.
610-1072	Sterile 15.2 x 244cm (6" x 96") accordion-folded cover with attached 3cm (1.2") latex-free NeoGuard tip with clips	12
610-833	Sterile 15.2 x 244cm (6" x 96") poly telescopically-folded cover with attached 3cm (1.2") latex-free NeoGuard tip with clips	12
610-834	Sterile 20.3 x 244cm (8" x 96") poly telescopically-folded cover with attached 6cm (2.4") latex-free NeoGuard tip with clips	12

Vist WWW.CIVCO.COM for a complete listing of surgical covers.



NEW!

Latex-Free Portable Ultrasound System Drape

- Provides a wide sterile field wherever ultrasound procedures are required
- Adhesive tape strips included with each system drape
- Safe means for reducing cross-contamination

REF. #	PRODUCT DESCRIPTION	QTY.
610-1037	Sterile 50.8 x 71.1cm (20" x 28") polyethylene system drape	24



SQUEEZE-MARK® Surgical Markers

- Unique squeezable barrel provides continuous flow for writing over wet prep solutions
- Available in highly visible, scrub-resistant, gentian violet ink

REF. #	PRODUCT DESCRIPTION	QTY.
610-746	Sterile SQUEEZE-MARK skin marking pen-medium tip	100
610-1108	Sterile SQUEEZE-MARK skin marking pen-medium tip (ruler not included)	100



Aquasonic® 100, Aquasonic Clear® & Scan®

- Aquasonic 100 is a high-viscosity coupling gel that is non-sensitizing, hypoallergenic and non-irritating
- Aquasonic Clear has the same formula and performance as Aquasonic 100, without perfume or color
- Scan is a less viscous gel and has a unique dual-viscosity formula which allows it to be applied as gel or fluid

REF. #	PRODUCT DESCRIPTION	QTY.
610-117	Aquasonic 100, 20 gram sterile gel packets (Ref. # 01-01)	48
610-130	Aquasonic 100, .25ml bottles (Ref. # 01-08)	12
610-116	Aquasonic 100, 5 liters, 1 dispenser bottle (Ref. # 01-50)	1
610-455	Aquasonic Clear, 5 liters econopack, 2 dispenser bottles (Ref. # 03-54)	4
610-892	Scan, 237ml bottle (Ref. # 11-08)	12
610-123	Scan, 1 gallon pump, 2 dispenser bottles (Ref. # 11-28S)	4



Ultra/Phonic™

- Long-lasting formula resists breakdown caused by salts in perspiration
- Will not run or dry out
- Available in sterile gel packets

REF. #	PRODUCT DESCRIPTION	QTY.
610-057	Ultra/Phonic, 20ml sterile gel packets (Ref. # NDC36-1000-20)	25



Cidex®, Cidex Plus®, Cidex OPA® & Enzol® Solutions

- Cidex activated dialdehyde solution is a fast and effective way to disinfect and sterilize many medical devices and instruments
- Cidex 14 day use destroys 99.8% of mycobacterium tuberculosis in 45 minutes at 25 degrees Celsius
- Cidex Plus (high-level disinfectant) is a 3.4% glutaraldehyde solution
- Cidex OPA (ortho-phthalaldehyde) is compatible with many medical devices, requires no activation or mixing, and destroys mycobacterium tuberculosis in 12 minutes at 20 degrees Celsius
- Enzol Enzymatic detergent removes tough, dried-on or hard-to-reach organic matter from instruments before disinfecting in Cidex

REF. #	PRODUCT DESCRIPTION	QTY.
610-640	Cidex, 14 day use, 1 gallon bottle (Ref. # 2250)	4
610-803	Cidex, 14 day use test strips, 15 strips per bottle (Ref. #2927)	2
610-641	Cidex Plus, 28 day use, 1 gallon bottle (Ref. # 2785)	4
610-804	Cidex Plus, 28 day use test strips, 15 strips per bottle (Ref. # 2926)	2
610-704	Cidex OPA, 14 day use, 1 gallon bottle (Ref. # 20390)	4
610-714	Cidex OPA test strips, 15 strips per bottle (Ref. # 20393)	2
610-643	Enzol, 1 gallon bottle (Ref. # 2252)	4



MetriCide® Sterilizing Solutions

- EPA registered as high-level disinfectant/sterilant solution
- High-level disinfection in 90 minutes at 25 degrees Celsius
- Simply activate and pour in container
- No dilution required
- Sterilization time is 10 hours at room temperature

REF. #	PRODUCT DESCRIPTION	QTY.
610-328	MetriCide 28, 28 day use, 1 quart bottle (Ref. # 10-2805)	16
610-805	MetriTest strips, 50 strips per bottle (Ref. # 10-304)	2



NEW!

Glutaraldehyde & OPA Spill Kits

- Contains Glute-Out® Neutralizer, an impervious gown, wrap around face shield, nitrile gloves, safety goggles and absorbent mats
- Large spill kits contain an added organic respirator, waterproof shoe covers plus additional Glute-Out Neutralizer and absorbent mats
- Meets OSHA and ANSI standards for safe use and handling of disinfectants in a healthcare facility
- Follow instructions for proper disposal of materials

REF. #	PRODUCT DESCRIPTION	QTY.
610-1098	Small spill kit with one .5oz bottle of Glute-Out to deactivate one quart of disinfectant (1), impervious gown (1), wrap-around face shield (1), nitrile gloves (1 pair), safety goggles (1) and large absorbent mats (2) (Ref. # AC-SP1)	1
610-1097	Large spill kit with one 2oz. bottle of Glute-Out to deactivate one gallon of disinfectant (1), impervious gown (1), wrap-around face shield (1), organic respirator (1), nitrile gloves (1 pair), safety goggles (1), waterproof shoe covers (1 pair), and large absorbent mats (4) (Ref. # AC-SP2)	1



NEW!

Glute-Out® Neutralizer for Glutaraldehyde & OPA

- Proven environmentally safe neutralizer
- Offers a solution to dangerous disinfectant disposal problems
- Neutralizes solutions in just five minutes
- Meets California EPA requirements

REF. #	PRODUCT DESCRIPTION	QTY.
610-1099	Glute-Out neutralizer, .5oz. bottles, each bottle deactivates one quart of disinfectant (Ref. # NG-QT)	24
610-1045	Glute-Out neutralizer, 2 oz. bottles (Ref. # NG-1GAL), each bottle deactivates one gallon of disinfectant	24

Visit WWW.CIVCO.COM for photos, pricing and system spec sheets.



Cleansing Towelettes

- Helps reduce risks associated with cross-contamination
- Pleasantly scented and treated with water, isopropyl alcohol and benzalkonium chloride
- Recommended for use as a general cleaner for hands and hard surfaces between procedures
- Contains 20% alcohol

REF. #	PRODUCT DESCRIPTION	QTY.
610-128	Towelettes, (7" x 6") 85 sheets per dispenser	6



Sani-Cloth HB® Germicidal Disposable Wipes

- Bactericidal, fungicidal and virucidal proven effective against 100 microorganisms
- Kills HBV, Influenza A-2 Virus, Acinetobacter, Aspergillus niger and Trichophyton mentagrophytes in 10 minutes
- Diamond embossed wipes contain less than 0.15% alcohol

REF. #	PRODUCT DESCRIPTION	QTY.
610-992	Sani-Cloth HB Germicidal Disposable Wipe, 20.3 x 35.6cm (8" x 14"), 65 sheets per dispenser (Ref. # Q85484)	6
610-991	Sani-Cloth HB Germicidal Disposable Wipe, 15.3 x 17.1cm (6" x 6.75"), 160 sheets per dispenser (Ref. # Q08472)	12



Sani-Cloth® Plus Germicidal Disposable Wipes

- Non-woven disposable cloths contain 14.85% alcohol
- Works in five minutes or less
- Proven to kill TB, MRSA, VRE, HIV-1, Staphylococcus aureus, Salmonella cholerasuis, Pseudomonas aeruginosa, Influenza A2 Hong Kong, Escherichia coli (pathogenic strain), and Herpes Simplex Type 2

REF. #	PRODUCT DESCRIPTION	QTY.
610-944	Sani-Cloth PLUS Germicidal Disposable Wipe, (8" x 14") 65 sheets per dispenser (Ref. # Q85084)	6
610-943	Sani-Cloth PLUS Germicidal Disposable Wipe, (6" x 6.75") 160 sheets per dispenser (Ref. # Q89072)	12



For more than 25 years CIVCO has been a trusted, single point of contact for ultrasound accessories

Complete Solutions Provider.

***CIVCO takes pride in customer satisfaction. As a valued CIVCO customer, we promise to provide you with high-quality products and unmatched personal service. We encourage you to let us know how we can better serve you.
Thank you for your business.***

Placing an Order

Corporate Office:
Toll-free within the U.S. and Canada:
(800) 445.6741
Toll-free Facsimile: (877) 329.2482
Direct Dial: (319) 656.4447
Facsimile: (319) 656.4451
E-mail: order@civco.com
WWW.CIVCO.COM

Office Hours

Monday-Friday, 7:00 a.m. - 5:30 p.m., CST.

To place an order after business hours, leave a voice mail telephone message or visit CIVCO's website. A customer service specialist will notify you on the next working day for order verification.

***Complete Warranty Available Online
at WWW.CIVCO.COM***

Payment Terms

Unless otherwise agreed by CIVCO in writing, all payments are due in full net thirty (30) days from date of invoice, in United States dollars, subject to approved credit at the time of order invoicing. Payment terms shall not be affected by any delay in installation or acceptance. Partial shipments will be billed as made and payable upon the terms above. No discount for early payment is authorized without CIVCO's written consent. CIVCO reserves the right to change its terms of sale or to require prior payment, letter of credit or COD when, in the opinion of CIVCO, the financial condition or previous payment record of any Purchaser so warrants. CIVCO shall not be obligated to extend credit to any Purchaser. Should a Purchaser become delinquent in any payment due, CIVCO shall be entitled, without notice, to suspend performance until the delinquency has been corrected. CIVCO shall retain a purchase money security interest in all Products sold, and proceeds thereof, until receipt of full payment in accordance with the terms hereof.

Transportation, Packing & Shipment Terms

Unless otherwise specified in CIVCO's order acknowledgement, Products will be shipped using CIVCO's designated carrier, CIP Destination, per Incoterms 2000, with all transportation, freight, insurance and similar charges, including all taxes, fees and import/export duties, at the expense of Purchaser. All such charges, taxes, and fees may be prepaid by CIVCO and added to Purchaser's invoice. Products shall be packaged for shipment in accordance with standard commercial practices; however, CIVCO shall not be responsible for delays or damages in shipment via a common carrier. Title to the Products and risk of loss shall pass to the Purchaser upon issuance of the carrier's bill of lading. CIVCO shall make reasonable efforts to deliver accepted orders promptly, and shall notify Purchasers of the reason for any delayed shipment. It is understood, however, that dates indicated for delivery or performance represent CIVCO's best

current estimates only and failure to perform on or by such dates shall not subject CIVCO to liability.

Quotations & Prices

All quotations furnished by CIVCO expire 90 days from the date thereof. All quotations constitute offers, and may be revoked or revised by CIVCO by written notice at any time prior to acceptance, except that quotations to non-U.S. customers are solicitations for offers to purchase. Clerical or typographical errors are subject to correction. Quoted prices and delivery dates are subject to product availability and receipt of prior orders, and prices are valid only for the quantities and payment schedule specified. Prices do not include, and Purchaser shall be responsible for, any applicable sales, use, excise, property, customs and similar taxes or duties. Purchaser shall also be responsible for itemized costs of transportation, special packaging, inspection charges, and insurance.

Claims, Shortages or Damages

Please open all shipments immediately and inspect for concealed damages, shortages or errors. If goods were damaged in shipping, please hold the package and request a carrier inspection and claim form. Please report any shortages or errors to your CIVCO customer service specialist within 30 days of receipt of shipment to arrange an authorized exchange.

Notice

Not all products available for international sale.

Sample Request

Free non-sterile samples are available on selected products. Samples are provided for sizing and fit purposes only and are not authorized for use in clinical procedures.

Return Policy

Product returns (excluding product recalls) shall be subject to CIVCO's product return policies, as in effect from time to time, and a minimum 20% restocking charge. All returned Product must be unused, in the original packaging, and in good condition with no defects or damage caused by Purchaser, its employees, officers, directors, agents, representatives or affiliates. Purchasers must obtain return authorization and instructions from CIVCO prior to returning any Product, including any Non-Conforming Product (any product not meeting CIVCO specifications as defined in CIVCO technical files). No Product may be returned if such Product has expired according to Product labels, if the seal or package integrity has been compromised, or if such Product or has not been used, handled or stored in accordance with CIVCO's quality control guidelines. Custom Product may be returned only if approved in writing by CIVCO and only upon payment of a reasonable restocking charge, which shall include but not be limited to tooling and development expenses incurred by CIVCO in connection with such Custom Product.

Privacy Policy

CIVCO's Commitment to Privacy Protection. Protecting the security and privacy of your personal information is important to CIVCO and to the way we conduct our business in compliance with laws on privacy, data protection and data security. We hope the policy outlined will help you understand what information CIVCO may collect, how CIVCO uses and safeguards that information and with whom we may share it.

When you provide CIVCO with personal information, it is used to respond to your inquiry, process your order or provide you access to specific account information. Also, is it used to support our customer relationship with you. The information is stored and processed and shared with CIVCO worldwide affiliates to better understand your business needs and how we can improve our products and services; or CIVCO may use personal information to contact you about a CIVCO offer in support of your business needs or to conduct online surveys to better understand our customers' needs. If you choose not to have your personal information used to support our customer relationship, we will respect your choice. We do not sell, rent or otherwise market your personal information to third parties.

Copyright Statement

All information provided in this catalog is protected by copyright. It is prohibited to copy, process, modify or commercially distribute this information without the express written permission of CIVCO.

The information herein provided is for general informative purposes only and no warranties or representations are made with respect thereto. Information in this catalog is prepared on the basis of the best information available at the time of publication. The catalog was accurate at the time of printing, but is subject to change as deemed appropriate by CIVCO. The information provided is subject to change at any time and without notice. Changed information may include, but is not limited to technical and commercial changes and changes in pricing. Binding declarations are only given after detailed inquiries.

Quality Standards

CIVCO achieved ISO certification in 1996. This international standard recognizes a commitment to quality in design, product development, manufacturing and customer service. CIVCO is currently a registered medical device establishment with the USA Food and Drug Administration and the CMDCAS-Canadian Medical Device Regulations. CIVCO is also certified to ISO 13485:2003 and Medical Device Directive 93/42 EEC requirements. Compliance with these regulations demonstrates CIVCO's commitment to providing the highest quality products and services to its customers. CIVCO's ISO certifications and assessments are performed on a regular basis by SGS United Kingdom Ltd. Systems & Services Certification.

COMPLETE SOLUTIONS PROVIDER

For more than 25 years, CIVCO has focused on designing, manufacturing and distributing products that improve patient care. Today our products help enhance life for patients undergoing radiation oncology, interventional radiology, ultrasound and urology services.

Contact us today to learn how you and your patients can benefit from a CIVCO solution.



For international distribution, visit CIVCO Worldwide at WWW.CIVCO.COM

Worldwide Headquarters: 102 First Street South, Kalona, Iowa 52247-9589 USA
Phone: 319.656.4447 | 800.445.6741 | Fax: 319.656.4451 | 877.329.2482
WWW.CIVCO.COM

COPYRIGHT © 2010. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ACCLISITE, ALIGN, CIVCO ASSIST, CIV-FLEX, COUPLED, CYBERMARK, ETRAX, HYDROSCAN, INFINITI, IVAS, MAGGI II PLUS, MULTI-PRO 2000, NEOGLARD, PRO-COVERS, QUICK-MOUNT, SURGI, SURGI-TIP, SURGI-BOOT, T-SHAPED SURGI, ULTRA-PRO II, ULTRA-PRO 3, ULTRA-PRO E, AND VIRTUTRAX ARE THE TRADEMARKS OF CIVCO. ALL OTHER TRADEMARKS ARE THE PROPERTY OF THEIR RESPECTIVE OWNERS. CIVCO PRODUCT NAMES MAY BE COVERED BY ONE OR MORE OF THE FOLLOWING U.S. PATENTS NOS. 5,941,889, 6,126,607, 6,361,499, 6,379,307 AND 7,395,563. OTHER PATENTS PENDING. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. PRINTED IN USA. 2010R-2466



Protex™ Disinfectant Spray

- Alcohol-free disinfectant
- Used on non-porous surfaces including transducers, mammography compressor plates and other non-surgical surfaces
- Provides protection against a broad spectrum of pathogens
- EPA approved and effective in killing a wide range of bacteria, viruses and fungi

REF. #	PRODUCT DESCRIPTION	QTY.
610-1079	Protex disinfectant spray, 12oz (Ref. # 42-12)	12



T-Spray™ & T-Spray II™

- Cleans and disinfects ultrasound probes and mammography compressor plates
- *T-Spray* will not cause transducer membrane swelling and will not leave grit or soap scum
- *T-Spray II* is an easy-to-use cleaner and disinfectant with no harsh solvents, alcohol, phenol or glutaraldehyde

REF. #	PRODUCT DESCRIPTION	QTY.
610-154	T-Spray transducer spray, 250ml pump spray (Ref. # NDC36-3200-25)	12
610-703	T-Spray II transducer spray, 250ml pump spray (Ref. # NDC36-2102-25)	12



DISPATCH® Disinfectant with Bleach

- Detergent and bleach solution (5500ppm NaOCl), stronger than 1:10 bleach solutions recommended by the Centers for Disease Control for disinfecting
- Kills staph, salmonella, Pseudomonas, human immunodeficiency virus, herpes simplex and polio virus in one minute
- Kills mycobacterium tuberculosis in two minutes

REF. #	PRODUCT DESCRIPTION	QTY.
610-705	DISPATCH liquid spray, 22 oz. (Ref. # 68967)	8



Transeptik®

- Economical, multi-purpose spray for use in cleaning ultrasound transducer surfaces between patient procedures
- Easy-to-use and will not cause transducer membrane swelling
- Transeptik does not contain chlorine, phenol or glutaraldehyde
- Accepted for use by most ultrasound equipment manufacturers

REF. #	PRODUCT DESCRIPTION	QTY.
610-388	Transeptik cleansing solution, 250ml pump spray (Ref. # 09-25)	12

REF. #	PAGE	REF. #	PAGE
610-001.....	6	610-932.....	10
610-002.....	6	610-943.....	14
610-004.....	6	610-944.....	14
610-005.....	6	610-946.....	10
610-013.....	6	610-989.....	6
610-015.....	6	610-991.....	14
610-020.....	7	610-992.....	14
610-022.....	7	610-1000.....	6
610-042.....	7	610-1010.....	6
610-043.....	7	610-1017.....	8
610-044.....	7	610-1018.....	8
610-046.....	7	610-1019.....	8
610-057.....	12	610-1030.....	9
610-116.....	12	610-1031.....	8
610-117.....	12	610-1032.....	8
610-123.....	12	610-1037.....	11
610-128.....	14	610-1045.....	13
610-130.....	12	610-1064.....	6
610-154.....	17	610-1072.....	11
610-208.....	7	610-1073.....	9
610-323.....	6	610-1074.....	9
610-328.....	13	610-1075.....	9
610-388.....	17	610-1079.....	17
610-455.....	12	610-1082.....	7
610-542.....	6	610-1083.....	7
610-575.....	6	610-1084.....	7
610-637.....	6	610-1088.....	7
610-640.....	12	610-1089.....	7
610-641.....	12	610-1090.....	7
610-643.....	12	610-1094.....	9
610-703.....	17	610-1095.....	9
610-704.....	12	610-1097.....	13
610-705.....	17	610-1098.....	13
610-714.....	12	610-1099.....	13
610-746.....	11	811-050.....	10
610-797.....	11	811-051.....	10
610-803.....	12	810-100.....	10
610-804.....	12	810-101.....	10
610-805.....	13	810-200.....	10
610-833.....	11	810-202.....	10
610-834.....	11	810-203.....	10
610-838.....	7	810-206.....	10
610-844.....	7	P09072.....	8
610-889.....	9	P09073.....	8
610-892.....	12	P09074.....	8



Want to place an order?

Call a Customer Service Specialist today for pricing or to place an order at 319.656.4447 or toll-free at 800.445.6741!