

PHILIPS C5-2 REFERENCE CHART

ULTRASOUND SYSTEM	TRANSDUCER	STARTER KIT	REPLACEMENT KIT
<ul style="list-style-type: none"> EnVisor - Only with Vision 2007 (C-I) upgrade 		 Ergo C5-2 Multi-Angle Ref.#: 667-071	 Guide: Ultra-Pro II TM Ref.#: 610-904 Gauge sizes: 8.5FR-23 (All sizes except 19)
<ul style="list-style-type: none"> EnVisor - Prior to Vision 2007 (C-I) upgrade HDI 4000 HDI 5000/5500cv HDI I iE33 iU22 		 Ergo C5-2 Single-Angle Ref.#: 613-241	 Guide: Multi-Pro 2000 TM Ref.#: 610-133 Gauge Sizes: 14, 16, 18, 20, 22, 25
<ul style="list-style-type: none"> HDI 1000/1500 HDI 3000/3500 		 C5-2 Single-Angle Ref.#: 2950-1517-02	 Guide: Multi-Pro 2000 TM Ref.#: 610-133 Gauge Sizes: 14, 16, 18, 20, 22, 25
<ul style="list-style-type: none"> HD3 		<div style="border: 1px solid blue; border-radius: 10px; padding: 10px; text-align: center;">Image not available</div> Ref.#: *BP-KIT-011-P (453561201041)	

* Not currently available through CIVCO. Contact Philips for availability.