

## SIEMENS EC9-4 & EV9-4 REFERENCE CHART

### ULTRASOUND SYSTEMS

### TRANSDUCER

### REUSABLE ENDOCAVITY NEEDLE GUIDE


### DISPOSABLE ENDOCAVITY NEEDLE GUIDES

### TRANSDUCERS COMPATIBLE WITH NEEDLE GUIDE

- G20
- G50
- G60 S
- X500(‡)
- Sequoia†



**REF. #:**  
• 610-666 (07848141)



**REF. #:**  
• 612-051  
• 612-061

- EC9-4
- EC7
- EC-10C5†
- EV9-4‡
- 6.5EC10

- G40
- X150
- X300



**REF. #:** 10041823



**REF. #:**  
• 676-132  
• 10041822

- EC9-4
- EV9-4

- G20
- G50
- G60 S



**REF. #:** 04392143



- EV9-4
- 6.5EV13
- 6.5EV13+
- 6.5EV13S

- Antares
- S2000



- EC9-4

\*Not currently available through CIVCO.  
Contact Siemens for product ordering information,  
888-826-9702.

**REF. #:**  
• 610-711  
• 610-842