


# Uni-frame® Thermoplastics

Uni-frame thermoplastics are one of the industry’s most popular choices. They provide the rigid fixation necessary for today’s radiotherapy treatment techniques. They are available with a disposable or reloadable frame and in a variety of perforation styles and thicknesses. A standard Uni-frame mask measures 9 x 10”; extended lengths are available for additional neck fixation. **Wide Uni-frame Masks** feature a thinner frame providing an expanded area of thermoplastic material to use.

**IMRT Reinforced Thermoplastics™** feature a combination of perforated and solid thermoplastic designed to provide rigidity. **ComfortPerf® thermoplastics** feature a special perforation pattern around the eyes and mouth, adding to patient comfort.



Patient in a standard perforated reloadable Uni-frame mask.




Disposable Thermoplastic (Left) vs. Reloadable Thermoplastic (Right)

Reloadable ABS Frame | MT-APUR-1




Reloadable frame with thermoplastic snapped out

## Uni-frame Mask Styles



Standard Perforated Uni-frame Mask



Wide Uni-frame Mask



ComfortPerf Uni-frame Mask



IMRT Reinforced Thermoplastic, Style 18, Uni-frame Mask



IMRT Reinforced Thermoplastic, Style 20, Uni-frame Mask



IMRT Reinforced Thermoplastic, Style 22, Uni-frame Mask

## Uni-frame® Thermoplas-

Perforation Type	Frame Type	Length	Thickness	Part Number
Standard Perforated	Disposable	10"	1.6mm	MT-APU-1.6
Standard Perforated	Disposable	10"	2.4mm	MT-APU-2.4
Standard Perforated	Disposable	10"	3.2mm	MT-APU
Standard Perforated	Reloadable	10"	1.6mm	MT-APUR-1.6
Standard Perforated	Reloadable	10"	2.4mm	MT-APUR-2.4
Standard Perforated	Reloadable	10"	3.2mm	MT-APUR
Standard Perforated	Disposable	12"	1.6mm	MT-APU-12-1.6
Standard Perforated	Disposable	12"	2.4mm	MT-APU-12-2.4
Standard Perforated	Disposable	12"	3.2mm	MT-APU-12
Standard Perforated	Reloadable	12"	1.6mm	MT-APUR-12-1.6
Standard Perforated	Reloadable	12"	2.4mm	MT-APUR-12-2.4
Standard Perforated	Reloadable	12"	3.2mm	MT-APUR-12
Standard Perforated	Disposable	14"	1.6mm	MT-APU-14-1.6
Standard Perforated	Disposable	14"	2.4mm	MT-APU-14-2.4
Standard Perforated	Disposable	14"	3.2mm	MT-APU-14
Standard Perforated	Reloadable	14"	1.6mm	MT-APUR-14-1.6
Standard Perforated	Reloadable	14"	2.4mm	MT-APUR-14-2.4
Standard Perforated	Reloadable	14"	3.2mm	MT-APUR-14
Standard Perforated	Disposable - Wide	10"	1.6mm	MT-APUW-1.6
Standard Perforated	Disposable - Wide	10"	2.4mm	MT-APUW-2.4
Standard Perforated	Disposable - Wide	10"	3.2mm	MT-APUW
Standard Perforated	Disposable - Wide	12"	2.4mm	MT-APUW-12-2.4
Standard Perforated	Disposable - Wide	12"	3.2mm	MT-APUW-12
ComfortPerf	Disposable	10"	2.4mm	MT-APU-CP201
ComfortPerf	Disposable	10"	3.2mm	MT-APU-CP301
ComfortPerf	Reloadable	10"	2.4mm	MT-APUR-CP201
ComfortPerf	Reloadable	10"	3.2mm	MT-APUR-CP301
IMRT Reinforced Style 18	Disposable	10"	2.4mm	MT-APUI-18-24
IMRT Reinforced Style 18	Disposable	10"	3.2mm	MT-APUI-18-32
IMRT Reinforced Style 18	Reloadable	10"	2.4mm	MT-APUIR-18-24
IMRT Reinforced Style 18	Reloadable	10"	3.2mm	MT-APUIR-18-32
IMRT Reinforced Style 20	Disposable	10"	2.4mm	MT-APUI-20-24
IMRT Reinforced Style 20	Disposable	10"	3.2mm	MT-APUI-20-32
IMRT Reinforced Style 20	Reloadable	10"	2.4mm	MT-APUIR-20-24
IMRT Reinforced Style 20	Reloadable	10"	3.2mm	MT-APUIR-20-32
IMRT Reinforced Style 22	Disposable	10"	3.2mm	MT-APUI-22-32
IMRT Reinforced Style 22	Reloadable	10"	3.2mm	MT-APUIR-22-32