

PCi⁺
MEDICAL
800-862-3394

Precautionary Statement
Dust may irritate eyes. Repeated exposure to skin may cause itching. If irritation occurs flush eyes with water and wash off skin in water or shower.
Read MSDS
Active Ingredient:
GLYCINE

Glute-Out[®] Neutralizer NG-1GAL

An environmentally friendly neutralizer for OPA & Glutaraldehyde

DOSE	INSTRUCTIONS FOR USE
1 bottle neutralizes 1 gallon (4 L) of OPA or Glutaraldehyde in 5 minutes.	Sprinkle into Glute/OPA. After 5 minutes dispose to drain while running cold water. Thoroughly rinse container.

NOTE: A dramatic color change will occur during activation.

Re-Order #
NG-1GAL

Northeast Silk Screen
O.Berk
PCi_NG-1Gal
5-3-12
10 x 16 305n