

# PEELSafe Advantage™ offers Convenience without Compromise

## Intuit Roll™ with PEELSafe Advantage



- Open face for easy insertion of transducer
- Pull-tab designed for consistent deployment without breaching sterile field
- Interior band remains intact, eliminating need for sliding additional bands down the cord

## Intuit SnapFit™ with PEELSafe Advantage



- Band snaps into place when cover is applied, eliminating the need for secondary band application
- Quick solution to replace tacky-residue wound dressing alternatives
- Open face for easy insertion of the transducer

REF #	PRODUCT DESCRIPTION	QTY.
610-2289	Intuit Roll sterile 14 x 147.3cm (5.5" x 58") rolled CIV-Flex cover with PEELSafe Advantage	24
610-2290	Intuit SnapFit sterile 14 x 30cm (5.5" x 12") accordion-folded CIV-Flex cover on applicator with PEELSafe Advantage	24



## Introducing Intuit™ Cover Line with PEELSafe Advantage™

As an infection control leader, CIVCO is excited to introduce a new transducer cover technology to save time and streamline your ultrasound workflow. PEELSafe Advantage offers fast and convenient application - simply peel off the label and apply directly to ultrasound transducer without any gel needed inside the cover.

The Intuit cover series combined with PEELSafe Advantage offers simplified application techniques, creating an improved user experience across all sterile ultrasound procedures. CIVCO's Intuit Series is an intuitive, easy-to-use sterile barrier. PEELSafe Advantage is designed to provide distortion-free images while scanning. This combination provides the end user convenience without compromise.



# PEELSafe

Advantage™

- Saves time with fast application
- Eliminates need for gel inside cover
- Quick and safe residue-free removal from transducer

