

The Right Way to Store Probes & Stay Compliant

“Store the device in a manner that will protect from damage or contamination... Hanging vertically in a cabinet and storing in a clean environment.”

-The Joint Commission



Actual Size:

24" (W) x 10" (D) x 74" (H)
61 x 25.4 x 188cm

Weight: 100 lbs. | 45.4kg

Minimum Clearance Required:

30" (W) x 36" (D) x 84" (H)
76.2 x 91.4 x 213.4cm

Safe & Clean Storage For Disinfected Ultrasound Probes

CIVCO's Storage Systems provide a safe and clean storage environment for your disinfected transesophageal and general purpose/endocavity ultrasound probes.

A 2022 study uncovered a **very low risk of pathogenic bacterium growth** on ultrasound probes **30 days** after high-level disinfection when stored in the CIVCO Storage Cabinet.¹ The cabinet's special design features:

- Fan Draws Air Through a Replaceable True HEPA Filter for Continuous Filtered Air Flow Using Positive Pressure
- Indicator Light Signals when Filter Replacement is Required
- Specialty Shelves Protect Delicate Probes From Damage
- Vertical Suspension of Probes Facilitates Drying
- Lockable Door Provides Secure and Safe Storage

REF #	PRODUCT DESCRIPTION	QTY
610-1323	TEE Probe Cabinet; Holds up to 3 probes; 100-240V; Includes U.S. power cord	1
610-1327	General Purpose & Endocavity Probe Cabinet; Holds up to 6 probes; 100-240V; Includes U.S. power cord	1
610-1328	Storage Cabinet replacement HEPA filter	1
610-2419	TEE Probe Cabinet; Holds up to 3 probes; 100-240V; Cord sold separately	1
610-2420	General Purpose & Endocavity Probe Cabinet; Holds up to 6 probes; 100-240V; Cord sold separately	1



CIVCO.com | 800.445.6741 | 319.248.6757

1. JOHNSON, R., SHARP, L., KIRK, B., WINTER, C., NATIVIDAD, M., MARIBEL, MEDINA, S., TAITE, L., DEAN, N., NELSON, C. (NOVEMBER/DECEMBER 2020). TRANSESOPHAGEAL ECHOCARDIOGRAM (TEE) PROBE TESTING: AN ANALYSIS TO EVALUATE THE FORMATION OF PATHOGENIC BACTERIAL COLONIES ON TEE PROBES IN STORAGE. RADIOLOGY MANAGEMENT, 43, 6, 9-17. COPYRIGHT © 2022. ALL RIGHTS RESERVED. CIVCO IS A REGISTERED TRADEMARK OF CIVCO MEDICAL SOLUTIONS. ALL PRODUCTS MAY NOT BE LICENSED IN ACCORDANCE WITH CANADIAN LAW. PRINTED IN USA. 2017R-3007 REV 4

LEARN MORE

