



GUS[®]

DISINFECTION
SOAK STATIONS[™]



-DISINFECT PROBES AND SCOPES OF ALL SIZES

-COST-EFFECTIVE HIGH-LEVEL DISINFECTION

Works Safely with Ortho-Phthalaldehyde (OPA), Hydrogen Peroxide and Glutaraldehyde

-ELIMINATE DUCTWORK & PLUMBING

-ALLOW FOR SAFE WORK ENVIRONMENT

Helps Meet The Joint Commission and Occupational Safety and Health Administration (OSHA) Standards



// G14KAL FOR GENERAL PURPOSE
AND ENDOCAVITY PROBES



// G14TC-3 FOR TEE PROBES

GUS® SYSTEMS FOR ULTRASOUND

	REF# (110v)	REF # (220v)	# OF PROBES	DESCRIPTION	TRAYS/TUBES
Endocavity or General Purpose	610-921	610-2141	1	GUS Disinfection Station G10VPL, Wall or Counter	(2) 3" x 17" (7.6 x 43.2cm)
	610-720	610-2150	2	GUS Disinfection Station G14KAL, Wall or Counter	(3) 3" x 17" (7.6 x 43.2cm)
Transesophageal (TEE or TOE)	610-850	610-2159	2	GUS Disinfection Station G14TC-3, Wall Mount	(3 Curved Tubes) 1.5" x 64" (3.8 x 162.6cm)

GUS Disinfection Soak Stations are not available for sale in every country. Customers outside the United States can contact distorder@civco.com or your local distributor with questions.

GUS® SYSTEMS FOR ENDOSCOPY

	REF# (110v)	REF# (220v)	DESCRIPTION	TRAYS/TUBES
Submersible Scopes	610-2165	610-2166	GUS Disinfection Station for Small Scopes G17, Space Saver, 19.5"H (49.5 cm H), Countertop	(2) 23.5" x 8.5" x 5.5" (59.7 x 22.2 x 14 cm)
	610-2168		GUS Disinfection Station for Small Scopes G18-E, Elevated Height, 24.5"H (62.2 cm H), Countertop	(2) 23.5" x 8.5" x 5.5" (59.7 x 22.2 x 14 cm)
	610-2173		GUS Disinfection Station for Large Scopes G32-E, Elevated Height, 24.5"H (62.2 cm H), Countertop	(2) 23.5" x 8.5" x 5.5" (59.7 x 22.2 x 14 cm)
	610-2174		GUS Disinfection Station for Large Scopes G32-S, Space Saver, 16.75"H (42.5 cm H), Countertop	(2) 23.5" x 8.5" x 5.5" (59.7 x 22.2 x 14 cm)



SAFETY ACCESSORIES - PROTECT YOURSELF & STAY COMPLIANT

FILTERS

Replace GUS filters every 6 months to effectively trap vapors from all high-level disinfectants.

REF#	DESCRIPTION	QTY
610-1207	Replacement Filter for GUS G10 Stations	2
610-1208	Replacement Filter for GUS G14 Stations	2
610-2102	Replacement Filter for GUS G17 Stations	2
610-2110	Replacement Filter for GUS G18 Stations	2
610-2125	Replacement Filter for GUS G32 Stations	2

DISINFECTANT NEUTRALIZER

Using a glycine-based powder, Glute Out® safely deactivates Glutaraldehyde and OPA, while Oxid-Out deactivates Hydrogen Peroxide. Visit CIVCO.com to view more available configurations.

REF#	DESCRIPTION	QTY
610-1099	Glute-Out Glutaraldehyde and OPA Neutralizer for 1 Quart (0.95 liters)	24
610-1045	Glute-Out Glutaraldehyde and OPA Neutralizer for 1 Gallon (3.8 liters)	24
610-2192	Glute-Out Glutaraldehyde and OPA Neutralizer for 5 Gallon (18.9 liters)	12
610-2202	Oxid-Out Hydrogen Peroxide Neutralizer for 1 Gallon (3.8 liters)	12
610-2203	Oxid-Out Hydrogen Peroxide Neutralizer for 8 Gallon (30 liters)	12
610-2204	Oxid-Out Hydrogen Peroxide Neutralizer for 1 Quart (0.95 liters)	24

DISINFECTANT WARMER KITS

Safely warm high-level disinfectant to meet minimum required disinfectant temperature.

REF#	DESCRIPTION	QTY
610-1291	Disinfectant Warmer for GUS Endocavity Stations, 110/120V - 17" (43.2 cm) Tubes - with G10VPL and G14KAL units	1

SPILL KIT

Meet Occupational Safety & Health Administration (OSHA) & American National Standards Institute (ANSI) guidelines with the spill kit.

REF#	DESCRIPTION	QTY
610-1098	Glute Out® Neutralizer (1), Organic Respirator (1), Nitrile Gloves (1), Impervious Gown (1), Bouffant Cap (1), Shoe Covers (1), Safety Goggles (1), Absorbent Mats (4), Water-Resistant Bucket (1)	1
610-2041	Spill Kit for Peracetic Acid & Hydrogen Peroxide (1), Organic Respirator (1), Nitrile Gloves (1), Impervious Gown (1), Bouffant Cap (1), Shoe Covers (1), Safety Goggles (1), Absorbent Mats (4), Water-Resistant Bucket (1)	1

STORAGE, TRANSPORT & ACCESSORIES

REF#	DESCRIPTION	QTY
610-1323	TEE Probe Storage Cabinet for up to 3 probes ; 100-240V; Includes U.S. power cord	1
610-1327	General Purpose & Endocavity Probe Storage Cabinet for up to 6 probes; 100-240V; Includes U.S. power cord	1
610-1418	TEE Probe Transport Tray	30
610-2037	Storage Rack for (3) Small Scopes	1
610-2016	Tube Cleaning Brush for G14TC-3 TEE Station	1
610-2049	Count Down Timer	1
610-2419	TEE Probe Storage Cabinet for up to 3 probes; 100-240V; Cordless - Not available in U.S. & Canada	1
610-2420	General Purpose & Endocavity Probe Storage Cabinet for up to 6 probes; 100-240V; Cordless - Not available in U.S. & Canada	1

Visit CIVCO.com or call **1-800-445-6741**

to learn more about GUS Disinfection Soak Stations and Safety Accessories