

Automate Your Ultrasound Probe Disinfection



Using ASTRA with Revital-Ox® RESERT® High Level Disinfectant

Ultrasound Probe Disinfection Guidelines

ANSI/AAMI ST91: 2015 Flexible and semi-rigid endoscope processing in health care facilities Section 5.7.1 General Considerations:

“High-level disinfection can be used in the reprocessing of semi-critical devices. Semi-critical devices can contact mucous membranes or non-intact skin and should preferably be sterilized or high-level disinfected with an FDA-cleared LCS/HLD prior to use on the next patient.”

Revital-Ox RESERT High Level Disinfectant

- Revital-Ox RESERT High Level Disinfectant is a **ready to use, fast, safe** hydrogen peroxide solution
- Up to 21 days reuse*
- **No additional ventilation** needed
- Leaves **no toxic residues** on devices
- No patient contraindications, contact dermatitis, or respiratory distress¹
- Biodegradable
- Safe to pour down drain¹
- No activation required



Revital-Ox® RESERT® R60 Solution Test Strip

- Dedicated test strip to test the minimum recommended concentration (MRC) of Revital-Ox RESERT High-Level Disinfectant
- Allows user to verify that the minimum recommended concentration (MRC) of hydrogen peroxide is present prior to each use
- Easy interpretation of results



ASTRA® Automated Reprocessor

for Transesophageal, Vaginal and Rectal Probe Disinfection

- **Automated data logging** and **consumable tracking** help you stay compliant with The Joint Commission standards. Clear, easy-to-follow steps guide you through the disinfection process
- Disinfect up to **two probes at a time** for faster turnaround and increased efficiency. 13- to 16-minute** total reprocessing time without sleep cycles
- Helps **reduce long-term operating costs** with reusable, off-the-shelf disinfectant compatibility



Product Descriptions

Revital-Ox RESERT HLD, 4X4L Bottle per case

Revital-Ox RESERT R60 Test Strip, 2 bottles, 60 test strips per bottle

ASTRA TEE Transesophageal Automated Probe Disinfection System

ASTRA VR Endovaginal/Endorectal Automated Probe Disinfection System



CIVCO.com/ASTRA | 1-800-445-6741

Contact CIVCO to learn more about disinfecting with ASTRA and Revital-Ox RESERT HLD and Revital-Ox RESERT R60 Solution Test Strip

*Conditional upon MRC passing result | **total time depends on ASTRA model used

References: 1. Revital-Ox® RESERT™ High Level Disinfectant; SDS No. 4455US [Online]; STERIS: Mentor, OH, January 18, 2018. <https://www1.steris.com/onbDocs/V434/1/970755.pdf>. (accessed June 9, 2020).