



# GUS<sup>®</sup>

DISINFECTION  
SOAK STATIONS<sup>™</sup>



**-DISINFECT PROBES AND SCOPES OF ALL SIZES**

**-COST-EFFECTIVE HIGH-LEVEL DISINFECTION**

*Works Safely with Ortho-Phthalaldehyde (OPA), Hydrogen Peroxide and Glutaraldehyde*

**-ELIMINATE DUCTWORK & PLUMBING**

**-ALLOW FOR SAFE WORK ENVIRONMENT**

*Helps Meet The Joint Commission and Occupational Safety and Health Administration (OSHA) Standards*



// G10VPL FOR ENDOCAVITY PROBES



// G14TC-3 FOR TEE PROBES

GUS® SYSTEMS FOR ULTRASOUND						
Endocavity or General Purpose	REF# (110v)	REF # (220v)	# OF PROBES	DESCRIPTION	TRAYS/TUBES	
	610-921	610-2141	1	GUS Disinfection Station G10VPL, Wall or Counter	(2) 3" x 17" (7.6 x 43.2cm)	
	610-720	610-2150	2	GUS Disinfection Station G14KAL, Wall or Counter	(3) 3" x 17" (7.6 x 43.2cm)	
Transesophageal (TEE or TOE)	610-850	610-2159	2	GUS Disinfection Station G14TC-3, Wall Mount	(3 Curved Tubes) 1.5" x 64" (3.8 x 162.6cm)	

GUS Disinfection Soak Stations are not available for sale in every country. Customers outside the United States can contact [distorder@civco.com](mailto:distorder@civco.com) or your local distributor with questions.

## GUS® SYSTEMS FOR ENDOSCOPY

	REF# (110v)	REF# (220v)	DESCRIPTION	TRAYS/TUBES
Non-Submersible Scopes	610-2144	610-2145	GUS Disinfection Station for 2 Small Scopes G14ENT-3, Wall Mount	(3) 1" x 19" (2.5 x 48.3 cm)
	610-2163		GUS Disinfection Station for 2 Urology Scopes G14-U3, Wall Mount	(3) 2" x 24" (5.1 x 61 cm)
Submersible Scopes	610-2165	610-2166	GUS Disinfection Station for Small Scopes G17, Space Saver, 19.5"H (49.5 cm H), Countertop	(2) 23.5" x 8.5" x 5.5" (59.7 x 22.2 x 14 cm)
	610-2168		GUS Disinfection Station for Small Scopes G18-E, Elevated Height, 24.5"H (62.2 cm H), Countertop	(2) 23.5" x 8.5" x 5.5" (59.7 x 22.2 x 14 cm)
	610-2173		GUS Disinfection Station for Large Scopes G32-E, Elevated Height, 24.5"H (62.2 cm H), Countertop	(2) 23.5" x 8.5" x 5.5" (59.7 x 22.2 x 14 cm)
	610-2174		GUS Disinfection Station for Large Scopes G32-S, Space Saver, 16.75"H (42.5 cm H), Countertop	(2) 23.5" x 8.5" x 5.5" (59.7 x 22.2 x 14 cm)
	610-2171		GUS Fume Hood G25-OER for Olympus Endoscope Reprocessors	



// G14ENT-3 FOR SMALL SCOPES



// G32-E FOR LARGE SCOPES

# SAFETY ACCESSORIES - PROTECT YOURSELF & STAY COMPLIANT

## FILTERS

Replace GUS filters every 6 months to effectively trap vapors from all high-level disinfectants.

REF#	DESCRIPTION	QTY
610-1207	Replacement Filter for GUS G10 Stations	2
610-1208	Replacement Filter for GUS G14 Stations	2
610-2102	Replacement Filter for GUS G17 Stations	2
610-2110	Replacement Filter for GUS G18 Stations	2
610-2119	Replacement Filter for GUS G25 Stations	2
610-2125	Replacement Filter for GUS G32 Stations	2

## DISINFECTANT NEUTRALIZER

Using a glycine-based powder, Glute Out® safely deactivates Glutaraldehyde and OPA, while Oxid-Out deactivates Peracetic Acid and Hydrogen Peroxide. Visit [CIVCO.com](http://CIVCO.com) to view more available configurations.

REF#	DESCRIPTION	QTY
610-1099	Glute-Out Glutaraldehyde and OPA Neutralizer for 1 Quart (0.95 liters)	24
610-2199	Glute-Out Caddy: Contains Wall Caddy, Glute-Out & Timer	1
610-2204	Oxid-Out Hydrogen Peroxide & Peracetic Acid Neutralizer for 1 Quart (0.95 liters)	24

## DISINFECTANT WARMER KITS

Safely warm high-level disinfectant to meet minimum required disinfectant temperature.

REF#	DESCRIPTION	QTY
610-1291	Disinfectant Warmer for GUS Endocavity Stations, 110/120V - 17" (43.2 cm) Tubes - with G10VPL and G14KAL units	1

## SPILL KIT

Spilled high-level disinfectant is dangerous. Meet Occupational Safety and Health Administration (OSHA) & American National Standards Institute (ANSI) guidelines with the spill kit.

REF#	DESCRIPTION	QTY
610-1098	Glute Out® Neutralizer (1), Organic Respirator (1), Nitrile Gloves (1), Impervious Gown (1), Bouffant Cap (1), Shoe Covers (1), Safety Goggles (1), Absorbent Mats (4), Water-Resistant Bucket (1)	1
610-2041	Spill Kit for Peracetic Acid & Hydrogen Peroxide (1), Organic Respirator (1), Nitrile Gloves (1), Impervious Gown (1), Bouffant Cap (1), Shoe Covers (1), Safety Goggles (1), Absorbent Mats (4), Water-Resistant Bucket (1)	1

## STORAGE, TRANSPORT & ACCESSORIES

REF#	DESCRIPTION	QTY
610-1323	TEE Probe Storage Cabinet for up to 3 probes	1
610-1327	General Purpose & Endocavity Probe Storage Cabinet for up to 6 probes	1
610-1418	TEE Probe Transport Tray	30
610-2037	Storage Rack for (3) Small Scopes	1
610-2016	Tube Cleaning Brush for G14TC-3 TEE Station	1
610-2020	Tube Cleaning Brush for G14ENT-3 Endoscope Stations	1
610-2049	Count Down Timer	1

Visit [CIVCO.com](http://CIVCO.com) or call **1-800-445-6741**  
to learn more about GUS Disinfection Soak Stations and Safety Accessories