

ULTRASOUND PROBES

Endovaginal & Endorectal Probes



1 Endocavity Probe
Wall or Counter
 17" Tubes (43.2 cm)

REF# 610-921: 110V Model
 REF# 610-2141: 220V Model



2 Endocavity Probes*
Wall or Counter
 17" Tubes (43.2 cm)

REF# 610-720: 110V Model
 REF# 610-2150: 220V Model

**Unit can be modified with REF# 610-2027 for use with general purpose probes*

Transesophageal Probes




2 Transesophageal Probes
Wall Mount
 1.5" x 64" Tubes (4 x 163 cm)

REF# 610-850: 110V Model
 REF# 610-2159: 220V Model

NON-SUBMERSIBLE ENDOSCOPES


ENT Scopes



2 ENT Scopes
Wall Mounted
 1" x 19" Tubes (2.5 x 48.3 cm)

REF# 610-2144: 110V Model
 REF# 610-2145: 220V Model

Cystoscopes



2 Cystoscopes
Wall Mounted
 2" x 24" Tubes (5.1 x 61 cm)

REF# 610-2163: 110V Model

GUS Disinfection Soak Stations are not available for sale in every country. Contact 1.800.445.6741 for more information.



SUBMERSIBLE ENDOSCOPES

ENT, Cystoscopes, and Hysteroscopes



Countertop

Space Saver 19.5"H (49.5 cm)
23.5"x 8.5"x 5.5" Trays (59.7 x 22.2 x 14 cm)
Stand-Alone Disposal System
Manual Two-Speed Fan

REF# 610-2165: 110V Model
REF# 610-2166: 220V Model



Countertop

Elevated Height 24.5"H (62.2 cm)
23.5" x 8.5" x 5.5" Trays (59.7 x 22.2 x 14 cm)
Built-In Disposal System
Automatic Two-Speed Fan

REF# 610-2168: 110V Model

Endoscopy, Multiple Endoscopes, Respiratory Therapy



Countertop

Elevated Height 24.5"H (62.2 cm)
15" x 21" x 6" Trays (38 x 53.3 x 15.2 cm) or
23.5" x 8.5" x 5.5" Trays (59.7 x 22.2 x 14 cm)
Built-In Disposal System
Automatic Two-Speed Fan

REF# 610-2173: 110V Model



Countertop

Space Saver 16.75"H (42.5 cm)
15" x 21" x 6" Trays (38 x 53.3 x 15.2 cm)
Built-In Disposal System
Automatic Two-Speed Fan

REF# 610-2174: 110V Model