

ULTRASOUND PROBES

Endovaginal & Endorectal Probes



1 Endocavity Probe
G10VPL Wall or Counter
 17" Tubes (43.2 cm)
 REF# 610-921: 110V Model
 REF# 610-2141: 220V Model



2 Endocavity Probes*
G14KAL Wall or Counter
 17" Tubes (43.2 cm)
 REF# 610-720: 110V Model
 REF# 610-2150: 220V Model
*Unit can be modified with REF# 610-2027 for use with general purpose probes

Transesophageal Probes



2 Transesophageal Probes
G14TC-3 Wall Mount
 1.5" x 64" Tubes (4 x 163 cm)
 REF# 610-850: 110V Model
 REF# 610-2159: 220V Model

ENDOSCOPES

ENT, Cystoscopes, and Hysteroscopes



G17 Countertop
 Space Saver 19.5"H (49.5 cm)
 23.5" x 8.5" x 5.5" Trays (59.7 x 22.2 x 14 cm)
 Stand-Alone Disposal System
 Manual Two-Speed Fan
 REF# 610-2165: 110V Model
 REF# 610-2166: 220V Model



G18-E Countertop
 Elevated Height 24.5"H (62.2 cm)
 23.5" x 8.5" x 5.5" Trays (59.7 x 22.2 x 14 cm)
 Built-In Disposal System,
 Automatic Two-Speed Fan
 REF# 610-2168: 110V Model

Endoscopy, Multiple Endoscopes, Respiratory Therapy



G32-E Countertop
 Elevated Height 24.5"H (62.2 cm)
 21" x 15" x 6" Trays (53.3 x 38 x 15.2 cm) or
 23.5" x 8.5" x 5.5" Trays (59.7 x 22.2 x 14 cm)
 Built-In Disposal System, Automatic Two-Speed Fan
 REF# 610-2173: 110V Model



G32-S Countertop
 Space Saver 16.75"H (42.5 cm)
 21" x 15" x 6" Trays (53.3 x 38 x 15.2 cm)
 Built-In Disposal System, Automatic Two-Speed Fan
 REF# 610-2174: 110V Model



GUS Disinfection Soak Stations are not available for sale in every country. Contact 1.800.445.6741 for more information.